



Fracture Table Foot Pad

- Use on leg-prep stand and fracture table foot holders.

40400	12 x 4 x ¼"	30.5 x 10 x .6 cm
--------------	-------------	-------------------



Knee Crutch Pads

- Three sizes for maximum versatility.

40402	Small	11 x 9 x ¾"	27.9 x 22.9 x .95 cm
40403	Medium	17 x 11 x ¼"	43 x 27.9 x .6 cm
40404	Large	22 x 20 x ½"	56 x 50 x 1.3 cm



Heel Protector Pad

- Fastens easily with convenient hook-and-loop attachments.

40401	12 ½ x 7 ¾ x ½"	32 x 20 x 1.3 cm
--------------	-----------------	------------------



Heel Support

- Secure and protect heel area without hook-and-loop attachment.
- Slit in back of cup cradles patient's Achilles tendon area.

40502	6 x 3 ¼" x 1 ¾"	15 x 8.3 x 4.5 cm
--------------	-----------------	-------------------



Heel Comfort Cup Set

- Use to complement anti-embolism stockings.
- Appropriate for recovery and radiology.
- Included AKTON strap fastens with convenient hook-and-loop attachments.

40520 (Set of 2)	4 ¾ x 4 x 1 ¾"	12.1 x 10.2 x 3.5 cm
-------------------------	----------------	----------------------



Resurfacing Kits

Quick and easy resurfacing of cuts or punctures on any ACTION Pad.

KIT002	4 x 4" 10 x 10 cm	KIT003	36 x 4" 91 x 10 cm
---------------	-------------------	---------------	--------------------

- VALUE product essential for any operating room.

