





- Action® Specialty Frame Pads provide comfort and stability for highly specialized surgeries.
- Shown in Photo: Wilson Frame, see below; OR Overlay Small, #40105; Supine Head Pad, #40200; Donut Head Pad, #40212; Chest rolls, #40620, #40621



<b>A</b>	<b>40701</b>	Wilson Frame Pad, 2-Piece Set (½" Thick)	29 x 10 x ½"	74 x 25.4 x 1.3 cm
<b>A</b>	<b>40702</b>	Wilson Frame Pad, 2-Piece Set (1" Thick)	29 x 10 x 1"	74 x 25.4 x 2.5 cm
<b>B</b>	<b>40705</b>	Andrews Frame™ Knee Pad	22 x 16 x 1"	56 x 40.6 x 2 cm
<b>B</b>	<b>40706</b>	Andrews Frame™ Iliac Strap	24 x 6 x ¾"	61 x 15 x .95 cm
<b>C</b>	<b>40110</b>	Montreal Frame Pad	18 ¼ x 12 x ½"	46.3 x 30.5 x 1.3 cm
<b>D</b>	<b>40712</b>	Kambin Spinal Frame Pad Set	34 x 14 ¼ x 8 ¾ x 10 ¾"	86.4 x 36.2 x 21.9 x 27.3 cm
<b>E</b>	<b>40711</b>	CF-94 Spinal Frame / CHOP	7 ½ diameter x 3" high	19 x 7.6 cm
<b>F</b>	<b>40700</b>	Relton Frame Pad Set, 4-Piece Set	8 ⅝ x 7 ⅝ x 4 ⅛" OD; 6 ½ x 5 ⅝ x 3 ⅝" ID	22.4 x 18.6 x 10.5 cm OD; 16.5 x 14.1 x 8.7 cm ID
<b>G</b>	<b>40703</b>	Fracture Table Hip Pad (Rectangular)	8 x 5 ½ x ½"	20 x 14 x 1.3 cm
<b>G</b>	<b>40707</b>	Fracture Table Hip Pad (Contoured)	1 ¾" diameter hole	4.5 cm diameter hole
<b>G</b>	<b>40704</b>	Fracture Table Post Pad	7 x 5 x ½"	18 x 12.7 x 1.3 cm
<b>H</b>	<b>40708</b>	Universal Lateral Positioner Iliac Crest Pad Set	3 x 3¼ x 1 ⅞" OD; 2 ½ x 2 ⅝ x 1" ID	7.6 x 8.3 x 4.8 cm OD; 5.5 x 6 x 2.5 cm ID
<b>H</b>	<b>40709</b>	Universal Lateral Positioner Frame Pad Cover Set	4 x 4 ½ x 2 ⅝" OD; 3 x 3 ¼ x 1 ⅞" ID	10 x 11.4 x 6 cm OD; 7.6 x 8.3 x 4.8 cm ID
<b>H</b>	<b>40710</b>	Universal Lateral Positioner Hip/Sacral Pad	26 x 10 x ½"	66 x 25.5 x 1.3 cm
<b>I</b>	<b>40100</b>	Action® Medium Overlay	46 x 20 x ½"	117 x 50 x 1.3 cm



"Positioning equipment should be used in a safe manner and according to manufacturer's written instructions."

2009 Perioperative Standards and Recommended Practices. RP: Positioning the Patient, AORN, Inc.

