

- Protects and cradles patient's head.
- Stabilizes head movement.
- Avoids neck overextension.



Donut Head Pads V

- Best seller.
- Variety of sizes for optimum positioning.
- Optional polymer lined dish provides premium protection.

40201	Adult	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40212	Adult Tall	8 OD x 3 ID x 3"	20.5 OD x 7.5 ID x 7.5 cm
40202	Pediatric	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40203	Neo-Natal	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40218	Adult Wide	8 OD x 4 1/2 ID x 1 3/4"	20.5 OD x 11.4 ID x 4.5 cm
40201D	Adult with Center Dish	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40202D	Pediatric with Center Dish	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40203D	Neo-Natal with Center Dish	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm



Horseshoe Head Pads

- For use in prone or lateral positions.

40204	Adult	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40211	Adult Tall	8 OD x 3 ID x 3"	20.5 OD x 7.5 ID x 7.5 cm
40205	Pediatric	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40206	Neo-Natal	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm

Ophthalmic Headrests

- Superior comfort and stability for procedures under local anesthesia.
- Surgeon arm/wrist support.

40209	Large	12 x 10 1/2 x 3 3/8"	31 x 27 x 10 cm
40210	Medium	12 x 10 x 2 1/8"	31 x 25 x 7.3 cm



Ophthalmic Cradle Headrests

- Reduced risk of alopecia.

40217	Adult	9 x 9 x 1 1/2"	23 x 23 x 3.8 cm
40217P	Pediatric	7 1/2 x 8 1/2 x 2 1/4"	19 x 21.6 x 5.7 cm



Lateral Head Pad with Center Dish

- Pre and Post-op use.
- Reduces pressure on side of patient's face and ear.



40207	9 x 10 x 1"	23 x 25 x 2.5 cm
--------------	-------------	------------------

Supine Head Pad

- AKTON polymer and foam construction.
- Built-in centering dish conforms for support and pressure prevention.



40200	9 x 10 x 2 1/4"	23 x 25 x 5.7 cm
--------------	-----------------	------------------

Prone Headrests

- Tube access via side or bottom channels.



40216M	Medium	8 3/4 x 7 3/4 x 4 3/4"	22.2 x 19.7 x 12 cm
40216L	Large	10 5/8 x 9 x 5 1/4"	27 x 22.0 x 13.4 cm

Contoured Head Pad

- Provides optimum stability of patient's head in the supine position.



40213	7 x 7 x 3 1/4" (1 1/2" contoured)	18 x 18 x 8 cm
--------------	-----------------------------------	----------------

Resurfacing Kits V

Quick and easy resurfacing of cuts or punctures on any ACTION Pad.

KIT002	4 x 4"	10 x 10 cm
KIT003	36 x 4"	91 x 10 cm

V - VALUE product essential for any operating room.

