

Pediatric Positioner Set

- Multipurpose set protects neonates to toddlers during a wide variety of procedures.
- Rectangular pad with indented channel allows chest expansion.
- Pediatric horseshoe pad (Model 40205) supports head in prone or lateral position.
- Includes two different size chest rolls (Model 40642 and 40646).

40114 Rectangular Pad: 20 x 9 x 2" 50.8 x 22.9 x 5 cm



Isolette Pad

- Consistent minimization of pressure and shear in NIC units to OR.

40120L Isolette Pad 26 x 16 x 1/2" 66 x 40.6 x 1.3 cm

Also available:

40120 Bassinet Pad 26 x 12 x 1/2" 66 x 30.5 x 1.3 cm



Pediatric Headrests

40202	Pediatric Donut Head Pad	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40203	Neo-Natal Donut Head Pad	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40202D	Pediatric Donut Head Pad with Center Dish	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40203D	Neo-Natal Donut Head Pad with Center Dish	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40205	Pediatric Horseshoe Head Pad	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40206	Neo-Natal Horseshoe Head Pad	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40217P	Ophthalmic Cradle Headrests Pediatric	7 1/2 x 8 1/2 x 2 1/4"	19 x 21.6 x 5.7 cm



Forehead Protector

- Protects infant, child, or adult patient from pressure exerted by oxygen tubes, etc.
- Semi-circular with flat bottom.

40208 2 1/2 x 1 3/8 x 3/8" 6.3 x 4 x 1.6 cm



Pediatric Chest Rolls

40639	6 x 2"	15 x 5 cm
40640	8 x 2"	20 x 5 cm
40641	10 x 2"	25.5 x 5 cm
40642	12 x 2"	30.5 x 5 cm
40643	6 x 3"	15 x 7.5 cm
40644	8 x 3"	20 x 7.5 cm
40645	10 x 3"	25.5 x 7.5 cm
40646	12 x 3"	30.5 x 7.5 cm



Resurfacing Kits

Quick and easy resurfacing of cuts or punctures on any ACTION Pad.

KIT002 4 x 4" 10 x 10 cm **KIT003** 36 x 4" 91 x 10 cm

- VALUE product essential for any operating room.

