



Ankle Strap Stirrup Pad Sets

- Adjustable to ensure perfect fit.
- AKTON polymer pads for lithotomy holders.
- Allows full view of feet.
- Includes all straps and hooks.



40510	Small	12 1/2 x 1 7/8 x 3/8"	31.7 x 4.7 x .95 cm
40511	Large	17 3/8 x 1 7/8 x 3/8"	44.1 x 4.7 x .95 cm



Boot/Stirrup Pad Set

- Protects lower leg and heel in Universal stirrups.

40500	10 1/2 x 6 x 14"	26.7 x 15.2 x 35.6 cm
--------------	------------------	-----------------------

Candy Cane Stirrup Pad Set

- Minimize risk of skin breakdown due to metal contact.
- Protects the peroneal nerve during lithotomy procedures.
- Eliminate tedious wrapping, unwrapping and re-wrapping.
- Use with Action Heel Protectors, Model 40401.

40501 (Set of 2)	28 x 2"OD	71 x 5 cm OD
-------------------------	-----------	--------------



SAM Boot Stirrup Pad Set

- AKTON polymer pockets protect even in radical elevated lithotomy.

40503	27 x 14 1/4T/5B x 1/4"	68.5 x 36T/12.5B x .6 cm
--------------	------------------------	--------------------------



Heel Protector Pad

- Fastens easily with convenient hook-and-loop attachments.

40401	12 1/2 x 7 3/4 x 1/2"	32 x 20 x 1.3 cm
--------------	-----------------------	------------------

- VALUE product essential for any operating room.

