

Maximum Decubitus and Tissue Trauma Protection in your O.R.
Clinically proven Action® O.R. Pads & Positioners

Action Products, Inc.

ACTION®

Body Smart® Solutions



www.actionproducts.com

Pressure Injury *Your Surgical Patients Are At Risk*

Overview

Pressure ulcer management has become a standard part of every modern hospital's protocol. The inclusion of HAPU's (hospital acquired pressure ulcers) as a Never Event by CMS in 2008 marked a turning point for most facilities.


As the pioneer of the polymer niche of pressure ulcer management, Action has been at the forefront of educating medical professionals and industry for over forty years. Our proprietary Akton® viscoelastic polymer still delivers the highest quality performance in hundreds of Body Smart® designs that feature its unique Shear Smart® surface. Designed to prevent capillary closure through the redistribution of pressure, the mitigation of shear forces and the even radiation of heat and humidity, our polymer significantly reduces the risk of tissue trauma. Few competitive products can claim the same.

Learn more about how to achieve your zero tolerance goals on our website featuring independent clinical studies, educational tools and programs for your staff at www.actionproducts.com.


How to Use This Catalog

Use this comprehensive catalog to source your facility's pressure management solutions. We have recently revised our catalog to reflect the critical and renewed emphasis on cost-containment within the hospital. We all know it is challenging to produce more with less. However, by moving to a reusable program of key, value-oriented products, your O.R. can do just that!



We have marked all our VALUE items with a  in an effort to draw your attention to the critical products key to providing a "basic" standard of care that complies with all the Patient Positioning Guidelines of the AORN. For your


convenience, we have compiled the entire list in the back of the catalog. Additionally, we offer a FREE Financial Comparative Study that proves the VALUE of savings achieved through transitioning to a reusable standard of care.



Look for this symbol at the bottom of each page identifying the products associated with preventing pressure ulcers in high risk areas.

Benefits of Action® Brand Products:

- Cost Effective
- Clinically Proven
- Will Not Support Bacterial Growth
- Biocompatible
- Imaging Compatible
- Fast and Easy Cleaning
- Latex Free
- Does Not Leak or Flow
- Tissue Equivalent
- Largest Product Selection
- Custom Manufacturing Available



Four major factors contributing to pressure ulcers:

- Uneven weight distribution
- Pressure
- Shear
- Heat and humidity build up



Hourglass Hand Table Overlay

- Appropriate for lengthy micro-hand procedures.

40330 32 x 24 x ½" 81.3 x 61 x 1.2 cm



Mayo Stand Cover

- Use for orthopedic procedures.
- Cushion costly instruments.
- Pillowcase style with rubber bottom.

40717 21 x 18 ½ x ⅝" 53.3 x 47 x 1.6 cm

Pillowcase



Fitted



Action® Disposable Overlay Covers

- Infection control product.
- Fitted sheet and pillowcase styles.
- Extend the service life of the pad.
- Improves room turnover rate.
- Case quantities available.

Pillowcase Style

40102 Made for table pad models 40100 & 40101 Case of 25

Fitted Sheet Style

40800 Made for table pad model 40100 Case of 25

40801 Made for table pad model 40101 Case of 25

40803 Made for table pad model 40103 Case of 25

40804 Made for table pad model 40104 Set of 3

40805 Made for table pad model 40105 Each

- Cited by AORN best practices for Pressure Ulcer Prevention.
- Provide pressure redistribution and reduce shear effects across entire table surface.
- Minimize patient movement.
- Low profile, simple design.
- ½" Akton® polymer for maximum effectiveness.



Action® O.R. Overlay with Perineal Cutout

- Ideal for lithotomy position.
- Wide range of applications in both regular and cysto operating rooms.

Action® O.R. Overlay, 3-Piece Set

- For use with segmented tables.
- Includes head, torso, and foot section pads.
- Provides pressure and shear reduction across entire table surface.

Action® O.R. Overlays

- Standard sizes and customs available.

40101	Large (Full-length)	72 x 20 x ½"	183 x 50 x 1.3 cm
40100	Medium	46 x 20 x ½"	117 x 50 x 1.3 cm
40105	Small - Hip pad	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm
40103	Medium with cutout	46 x 20 x ½" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
40104A	Head	9 x 20 x ½"	22.9 x 50.8 x 1.3 cm
40104B	Torso	39 ½ x 20 x ½"	100.3 x 50.8 x 1.3 cm w/ perineal cutout
40104C	Foot	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm w/ perineal cutout



Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.


KIT002 4 x 4" 10 x 10 cm **KIT003** 36 x 4" 91 x 10 cm

☑ - VALUE product essential for any operating room.




ACTION® PERFORMANCE SERIES O.R. TABLE PADS

- 4-way stretch antimicrobial fabric surface.
- Continuous layer of 3/8" Akton® polymer.
- Soft viscofoam mid-layer.
- High-density foam base.
- Fluid proof, welded seams.
- 100% no-stitch hook-and-loop attachment.
- Proprietary hidden vent system.
- Easy to clean waterfall construction.
- Available in multiple thicknesses.
- Armboards available.
- Custom pads are available. Call for details.



3" Performance Pad

5'11" Male Subject, 175 lbs. Average Pressure (mmHg): 18.72 Room Temp: 68 degrees Fahrenheit
 Equipment: FSA Pressure Map System Range: 1 – 150 mmHg



Introducing our **Performance** and **Performance PLUS Series** of table pads that will help you achieve your zero tolerance goals for HAPUs (Hospital Acquired Pressure Ulcers) and fall well within your budget.

Our pads are clinically engineered to give the maximum in decubiti prevention, infection control, and stability for the ultimate in patient safety. With two styles and multiple thicknesses, we have the perfect mattress for any budgetary requirement.

Action® Performance Series O.R. Table Pads


55120	2" Fabric Top	Fits Amsco 3080/3085
55130	2" Fabric Top	Fits Steris C-Max
55310	2" Fabric Top	Fits Skytron 3500
55320	2" Fabric Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
55330	2" Fabric Top	Fits Skytron 6700/ 6701
55370	2" Fabric Top	Fits Skytron 3501
57120	3" Fabric Top	Fits Amsco 3080/3085
57130	3" Fabric Top	Fits Steris C-Max
57310	3" Fabric Top	Fits Skytron 3500
57320	3" Fabric Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
57330	3" Fabric Top	Fits Skytron 6700/ 6701
57370	3" Fabric Top	Fits Skytron 3501
59100	2" Fabric Top	Universal Armboard
59600	3" Fabric Top	Universal Armboard




Durable, welded, midline seam and proprietary hidden vent structure enhance infection control for maximum pad life.



ACTION® PERFORMANCE PLUS SERIES O.R. TABLE PADS



- Shear Smart® surface of exposed 3/8" Akton® polymer.
- Multi-directional stretch film covering on surface.
- Soft viscofoam mid-layer.
- High-density foam base.
- Tissue equivalent.
- Does not support bacterial growth.
- Fluid proof, welded seams.
- 100% no-stitch hook-and-loop attachment.
- Proprietary hidden vent system.
- Easy to clean waterfall construction.
- Available in multiple thicknesses.
- Armboards available.
- Custom pads are available. Call for details.



3" Performance Plus Pad

5' 11" Male Subject, 175 lbs. Average Pressure (mmHg): 22.33 Room Temp: 69 degrees Fahrenheit
Equipment: FSA Pressure Map System Range: 1 – 150 mmHg

Action® Performance Plus Series O.R. Table Pads

50120	2" Polymer Top	Fits Amsco 3080/3085
50130	2" Polymer Top	Fits Steris C-Max
50310	2" Polymer Top	Fits Skytron 3500
50320	2" Polymer Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
50330	2" Polymer Top	Fits Skytron 6700/ 6701
50370	2" Polymer Top	Fits Skytron 3501
52120	3" Polymer Top	Fits Amsco 3080/3085
52130	3" Polymer Top	Fits Steris C-Max
52310	3" Polymer Top	Fits Skytron 3500
52320	3" Polymer Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
52330	3" Polymer Top	Fits Skytron 6700/ 6701
52370	3" Polymer Top	Fits Skytron 3501
54100	2" Polymer Top	Universal Armboard
54600	3" Polymer Top	Universal Armboard



Akton® polymer in a proprietary soft film that reduces shear effects for optimal patient safety.



Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

KIT002 4 x 4" 10 x 10 cm

KIT003 36 x 4" 91 x 10 cm

- VALUE product essential for any operating room.




HEAD PADS

- Protects and cradles patient's head.
- Stabilizes head movement.
- Avoids neck overextension.

Lateral Head Pad with Center Dish

- Pre and Post-op use.
- Reduces pressure on side of patient's face and ear.



40207 9 x 10 x 1" 23 x 25 x 2.5 cm

Supine Head Pad

- Akton® polymer and foam construction.
- Built-in centering dish conforms for support and pressure prevention.



40200 9 x 10 x 2 1/4" 23 x 25 x 5.7 cm

Prone Headrests

- Tube access via side or bottom channels.



40216M Medium 8 3/4 x 7 3/4 x 4 3/4" 22.2 x 19.7 x 12 cm

40216L Large 10 3/8 x 9 x 5 1/4" 27 x 22.0 x 13.4 cm

Contoured Head Pad

- Provides optimum stability of patient's head in the supine position.



40213 7 x 7 x 3 1/4" (1 1/2" contoured) 18 x 18 x 8 cm

Resurfacing Kits V

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

KIT002 4 x 4" 10 x 10 cm

KIT003 36 x 4" 91 x 10 cm



V - VALUE product essential for any operating room.

Donut Head Pads V

- Best seller.
- Variety of sizes for optimum positioning.
- Optional polymer lined dish provides premium protection.



40201 Adult 8 OD x 3 ID x 1 1/4" 20.5 OD x 7.5 ID x 4.5 cm

40212 Adult Tall 8 OD x 3 ID x 3" 20.5 OD x 7.5 ID x 7.5 cm

40202 Pediatric 5 1/2 OD x 2 1/4 ID x 1 1/4" 14 OD x 5.7 ID x 3.2 cm

40203 Neo-Natal 3 1/4 OD x 1 1/2 ID x 3/4" 8.3 OD x 3.8 ID x 2 cm

40218 Adult Wide 8 OD x 4 1/2 ID x 1 3/4" 20.5 OD x 11.4 ID x 4.5 cm

40201D Adult with Center Dish 8 OD x 3 ID x 1 1/4" 20.5 OD x 7.5 ID x 4.5 cm

40202D Pediatric with Center Dish 5 1/2 OD x 2 1/4 ID x 1 1/4" 14 OD x 5.7 ID x 3.2 cm

40203D Neo-Natal with Center Dish 3 1/4 OD x 1 1/2 ID x 3/4" 8.3 OD x 3.8 ID x 2 cm

Horseshoe Head Pads

- For use in prone or lateral positions.



40204 Adult 8 OD x 3 ID x 1 1/4" 20.5 OD x 7.5 ID x 4.5 cm

40211 Adult Tall 8 OD x 3 ID x 3" 20.5 OD x 7.5 ID x 7.5 cm

40205 Pediatric 5 1/2 OD x 2 1/4 ID x 1 1/4" 14 OD x 5.7 ID x 3.2 cm

40206 Neo-Natal 3 1/4 OD x 1 1/2 ID x 3/4" 8.3 OD x 3.8 ID x 2 cm

Ophthalmic Headrests

- Superior comfort and stability for procedures under local anesthesia.
- Surgeon arm/wrist support.



40209 Large 12 x 10 1/2 x 3 1/8" 31 x 27 x 10 cm

40210 Medium 12 x 10 x 2 1/8" 31 x 25 x 7.3 cm

Ophthalmic Cradle Headrests

- Reduced risk of alopecia.



40217 Adult 9 x 9 x 1 1/2" 23 x 23 x 3.8 cm

40217P Pediatric 7 1/2 x 8 1/2 x 2 1/4" 19 x 21.6 x 5.7 cm

ARM / SHOULDER / ELBOW PADS



Contoured Armboard Pads

- Full-arm protection: biceps, elbows, and forearms.
- Contour minimizes chance of rolling patient's arm.
- Ideal for prone positioning, especially laminectomy procedures.

40307	Small	16 x 5 1/2 x 2"	41 x 14 x 5 cm
40310	Medium	20 x 5 3/4 x 2 1/4"	51 x 15 x 6 cm
40306	Large	24 x 5 3/4 x 2 1/4"	61 x 15 x 6 cm



Armboard Pads V

- Protect patient's bicep, elbow, and forearm areas.
- Use on standard armboard table extensions.

40300	Short	15 x 4 x 1/4"	38 x 10 x .6 cm
40301	Long	24 x 4 1/2 x 3/8"	61 x 11.5 x .95 cm
40309	Wide Short	15 x 6 1/2 x 1/4"	38 x 16.5 x .6 cm
40308	Wide Long	24 x 6 1/2 x 3/8"	61 x 16.5 x .95 cm



Hand/Wrist Support

NEW all Akton® polymer construction.

- Superior hand and wrist stabilization and protection.
- Perfect for Carpal Tunnel, Arterial Cannulization, Dialysis, and other procedures.
- Ideally positions hand/wrist for harvesting Radial Artery during Coronary Artery Bypass surgery.
- Best used with stretch bandages for stability.

40321	11 x 3 1/4 x 2 3/8"	28 x 8.3 x 6 cm
--------------	---------------------	-----------------



Toboggan Pad

- Use pad in conjunction with arm/leg toboggan shield.
- Prevents patient's arm or leg from extending over side of table.
- Protects against leaning against patient's arm or leg.
- Reduces risk from shear or pressure from direct contact with shield.
- Set includes the shield.

40130	Toboggan Pad	16 1/4 x 8 1/4 x 3/8"	41.3 x 21 x .95 cm
40130-S*	Toboggan Pad Set	16 1/4 x 8 1/4 x 3/8"	41.3 x 21 x .95 cm

*Warranty covers pad only.



Arm/Hand Stabilizer Strap

- Stabilizes arm on the armboard.
- Fast, easy hook-and-loop fasteners.

40315	25 x 2 x 1/4"	63.5 x 5 x .635 cm
--------------	---------------	--------------------

Ulnar Forearm Protector

- Protects patient's forearm.
- Convenient, quick hook-and-loop attachments.

40303	18 x 13 x 1/4"	46 x 33 x .6 cm
--------------	----------------	-----------------



Ulnar Nerve Protector V

- Protects ulnar nerve while allowing full access to arm.
- Convenient, quick hook-and-loop attachments.

40304	Small	15 3/4 x 6 x 1/2"	40 x 15 x 1.3 cm
40304L	Large	19 x 8 1/2 x 1/2"	48 x 21.6 x 1.3 cm



Shoulder Brace, 2-Piece Set

- Protect shoulders in Trendelenburg position.
- Pads wrap around braces with hook-and-loop attachments.

40302	5 x 6 x 1"	13 x 15 x 2.5 cm
--------------	------------	------------------

V - VALUE product essential for any operating room.



ANKLE / FOOT / HEEL / KNEE PADS



Fracture Table Foot Pad

- Use on leg-prep stand and fracture table foot holders.

40400	12 x 4 x ¼"	30.5 x 10 x .6 cm
--------------	-------------	-------------------



Knee Crutch Pads V

- Three sizes for maximum versatility.

40402	Small	11 x 9 x ¾"	27.9 x 22.9 x .95 cm
40403	Medium	17 x 11 x ¼"	43 x 27.9 x .6 cm
40404	Large	22 x 20 x ½"	56 x 50 x 1.3 cm



Heel Comfort Cup Set

- Use to complement anti-embolism stockings.
- Appropriate for recovery and radiology.
- Included polymer strap fastens with convenient hook-and-loop attachments.

40520 (Set of 2)	4 ¾ x 4 x 1 ⅝"	12.1 x 10.2 x 3.5 cm
-------------------------	----------------	----------------------



Heel Support V

- Secure and protect heel area without hook-and-loop attachment.
- Channel in back of cup cradles patient's Achilles tendon area.

40502	6 x 3 ¼" x 1 ¾"	15 x 8.3 x 4.5 cm
--------------	-----------------	-------------------



Heel Protector Pad

- Fastens easily with convenient hook-and-loop attachments.

40401	12 ½ x 7 ¾ x ½"	32 x 20 x 1.3 cm
--------------	-----------------	------------------



V - VALUE product essential for any operating room.



Ankle Strap Stirrup Pad Sets

- Adjustable to ensure perfect fit.
- Akton® polymer pads for lithotomy holders.
- Allows full view of feet.
- Includes all straps and hooks.


40510	Small	12 ½ x 1 ½ x ¾"	31.7 x 4.7 x .95 cm
40511	Large	17 ¾ x 1 ½ x ¾"	44.1 x 4.7 x .95 cm



Candy Cane Stirrup Pad Set

- Minimize risk of skin breakdown due to metal cane contact.
- Protects the peroneal nerve during lithotomy procedures.
- Eliminates tedious wrapping, unwrapping and re-wrapping.
- Use with Action® Heel Protectors, Model 40401.

40501 (Set of 2)	28 x 2"OD	71 x 5 cm OD
-------------------------	-----------	--------------



Boot/Stirrup Pad Set

- Protects lower leg and heel in Universal stirrups.

40500	10 ½ x 6 x 14"	26.7 x 15.2 x 35.6 cm
--------------	----------------	-----------------------



SAM Boot Stirrup Pad Set

- AKTON polymer pockets protect even in radical elevated lithotomy.

40503	27 x 14 ¼T/5B x ¼"	68.5 x 36T/12.5B x .6 cm
--------------	--------------------	--------------------------

Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

KIT002	4 x 4"	10 x 10 cm
KIT003	36 x 4"	91 x 10 cm

- VALUE product essential for any operating room.



POSITIONING CHEST ROLLS



Contoured Chest Rolls

- Conforms to an adult's rib cage.

40622	Small	12 x 6 ¼ x 5"	30.5 x 16 x 13 cm
40623	Medium	16 x 6 ¼ x 5"	40.6 x 16 x 13 cm
40624	Large	20 x 6 ¼ x 5"	50.8 x 16 x 13 cm



Chest Positioners

- Conforms to patient's shoulder-to-chest area.
- Rounded edges for axilla.
- Helps alleviate pressure on thorax area.

40625	Left	28/23 x 8 x 3"	71/58.5 x 20.3 x 7.6 cm
40626	Right	28/23 x 8 x 3"	71/58.5 x 20.3 x 7.6 cm



Roll Covers V

- Wraps around stapl check-type rolls your facility already has.
- Upgrade basic chest rolls for pressure redistribution.
- Features hook-and-loop attachments.

40600	Chest Roll Cover	18 x 18 x ½"	46 x 46 x 1.3 cm
40601	Axillary Roll Cover	17 x 12 x ½"	43 x 30.5 x 1.3 cm

- Increased stability.
- Provides support for upper body positioning.



Flat-Bottomed Chest Rolls


- Stable, no-slip design.
- Protects entire chest area from pressure in prone position.
- Wide range of adult and pediatric sizes to accommodate virtually any need.

Adult Chest Rolls			Pediatric Chest Rolls		
40616	12 x 4"	30.5 x 10 cm	40639	6 x 2"	15 x 5 cm
40617	16 x 4"	40.6 x 10 cm	40640	8 x 2"	20 x 5 cm
40618	20 x 4"	50.8 x 10 cm	40641	10 x 2"	25.5 x 5 cm
40619	12 x 6"	30.5 x 15 cm	40642	12 x 2"	30.5 x 5 cm
40620	16 x 6"	40.6 x 15 cm	40643	6 x 3"	15 x 7.5 cm
40621	20 x 6"	50.8 x 15 cm	40644	8 x 3"	20 x 7.5 cm
			40645	10 x 3"	25.5 x 7.5 cm
			40646	12 x 3"	30.5 x 7.5 cm


Resurfacing Kits V

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

KIT002	4 x 4"	10 x 10 cm	KIT003	36 x 4"	91 x 10 cm
---------------	--------	------------	---------------	---------	------------



V - VALUE product essential for any operating room.



Trapezoid Positioner

- Provides clear leg separation.

40630 10 x 7 x 3½" 25.4 x 17.8 x 8.9 cm

Patient Positioners

- Versatile positioners eliminate use of old-fashioned sandbags.
- Polymer positioners are a good solution to unsafe rolled towels, IV bags and sandbags.

40602 Small 10 x 4 ½ x 2" 25.4 x 11.4 x 5 cm

40605 Large 14 x 4 ½ x 2" 35.5 x 11.4 x 5 cm


40615 Extra Long 18 x 4 x 2" 46 x 10 x 5 cm



Dome-Shaped Positioner

- Provides pressure protection under the knees in supine position.


40603 14 x 6 x 3" 36 x 15 x 7.6 cm



Lateral Positioner

- Accommodates both shoulder and upper arm.
- Allows for pressure redistribution and nerve protection.
- Lightweight foam base with AKTON polymer pad overlay.

40610 28 x 19.5 x 6/1" 71.1 x 49.5 x 15.2/2.5 cm



Action® Positioning Strap

- Use to minimize patient movement.
- Hook-and-loop attachments connect to rails on tables.
- Akton® polymer adhered to black elastomeric band.
- Akton® polymer area measures 38 ½ x 3 ¼ x ¼" (98 x 8.2 x .6 cm).

40713 56½ x 3½ x ¼" 143.5 x 8.9 x .6 cm

Graduated Thoracic Positioner

- Tapered design allows for full chest extension.
- Designed for cardiovascular or thoracic procedures.

40604 14 x 9 x 2/1 ½" 36 x 23 x 5/3.8 cm



Vac-Pac® Covers

- Protect patient tissue integrity by laying this pad on top of a Vac-Pac® prior to compression.
- Minimize shear and pressure created by the harsh vacuum positioner.
- Ideal for total hip and kidney/liver procedures.

40612 Universal Vac-Pac® Cover 37 x 28 x ½" 94 x 71 x 1.3 cm

40613 T-Shaped Vac-Pac® Cover 34 x 15 x ½" bottom 86.5 x 38 x 1.3 cm bottom
23 x 9 x ½" top 58.5 x 23 x 1.3 cm top

Vac-Pac is a registered trademark of Olympic Medical.

- VALUE product essential for any operating room.



PEDIATRIC POSITIONERS

Pediatric Positioner Set V

- Multipurpose set protects neonates to toddlers during a wide variety of procedures.
- Rectangular pad with indented channel allows chest expansion.
- Pediatric horseshoe pad (Model 40205) supports head in prone or lateral position.
- Includes two different size chest rolls (Model 40642 and 40646).

40614 Rectangular Pad: 20 x 9 x 2" 50.8 x 22.9 x 5 cm



Isolette Pad

- Consistent minimization of pressure and shear in NIC units to O.R.

40120L Isolette Pad 26 x 16 x 1/2" 66 x 40.6 x 1.3 cm

Also available:

40120 Bassinet Pad 26 x 12 x 1/2" 66 x 30.5 x 1.3 cm



Pediatric Headrests

40202	Pediatric Donut Head Pad	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40203	Neo-Natal Donut Head Pad	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40202D	Pediatric Donut Head Pad with Center Dish	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40203D	Neo-Natal Donut Head Pad with Center Dish	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40205	Pediatric Horseshoe Head Pad	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40206	Neo-Natal Horseshoe Head Pad	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40217P	Ophthalmic Cradle Headrests Pediatric	7 1/2 x 8 1/2 x 2 1/4"	19 x 21.6 x 5.7 cm

For more sizes please view page 6.



Forehead Protector

- Protects infant, child, or adult patient from pressure exerted by oxygen tubes, etc.
- Semi-circular with flat bottom.

40208 2 1/2 x 1 5/8 x 5/8" 6.3 x 4 x 1.6 cm



For more sizes please view page 10.

Pediatric Chest Rolls

40639	6 x 2"	15 x 5 cm
40640	8 x 2"	20 x 5 cm
40641	10 x 2"	25.5 x 5 cm
40642	12 x 2"	30.5 x 5 cm
40643	6 x 3"	15 x 7.5 cm
40644	8 x 3"	20 x 7.5 cm
40645	10 x 3"	25.5 x 7.5 cm
40646	12 x 3"	30.5 x 7.5 cm



V - VALUE product essential for any operating room.



Resurfacing Kits V

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.






KIT002 4 x 4" 10 x 10 cm **KIT003** 36 x 4" 91 x 10 cm

SPECIALTY FRAME PADS

- Action® Specialty Frame Pads provide comfort and stability for highly specialized surgeries.
- Shown in Photo: Wilson Frame, see below; O.R. Overlay Small, p.3; Supine Head Pad, p.6; Donut Head Pad, p.6; Chest Rolls, p.10



A	40701	Wilson Frame Pad, 2-Piece Set (½" Thick)	29 x 10 x ½"	74 x 25.4 x 1.3 cm
A	40702	Wilson Frame Pad, 2-Piece Set (1" Thick)	29 x 10 x 1"	74 x 25.4 x 2.5 cm
B	40705	Andrews Frame™ Knee Pad	22 x 16 x 1"	56 x 40.6 x 2 cm
B	40706	Andrews Frame™ Iliac Strap	24 x 6 x ¾"	61 x 15 x .95 cm
C	40110	Montreal Frame Pad	18 ¼ x 12 x ½"	46.3 x 30.5 x 1.3 cm
D	40712	Kambin Spinal Frame Pad Set	34 x 14 ¼ x 8 ⅝ x 10 ¾"	86.4 x 36.2 x 21.9 x 27.3 cm
E	40711	CF-94 Spinal Frame / CHOP	7 ½ diameter x 3 high x ⅝"	19 x 7.6 x 1.6 cm
F	40700	Relton Frame Pad Set, 4-Piece Set	8 ⅝ x 7 ⅝ x 4 ¼" OD; 6 ½ x 5 ⅝ x 3 ⅜" ID	22.4 x 18.6 x 10.5 cm OD; 16.5 x 14.1 x 8.7 cm ID
G	40703	Fracture Table Hip Pad (Rectangular)	8 x 5 ½ x ½"	20 x 14 x 1.3 cm
G	40707	Fracture Table Hip Pad (Contoured)	1 ¾" diameter hole	4.5 cm diameter hole
G	40704	Fracture Table Post Pad	7 x 5 x ½"	18 x 12.7 x 1.3 cm
H	40708	Universal Lateral Positioner Iliac Crest Pad Set	3 x 3 ¼ x 1 ⅞" OD; 2 ⅞ x 2 ⅜ x 1" ID	7.6 x 8.3 x 4.8 cm OD; 5.5 x 6 x 2.5 cm ID
H	40709	Universal Lateral Positioner Frame Pad Cover Set	4 x 4 ½ x 2 ⅝" OD; 3 x 3 ¼ x 1 ⅞" ID	10 x 11.4 x 6 cm OD; 7.6 x 8.3 x 4.8 cm ID
H	40710	Universal Lateral Positioner Hip/Sacral Pad	26 x 10 x ½"	66 x 25.5 x 1.3 cm
I	40100	Action® Medium Overlay	46 x 20 x ½"	117 x 50 x 1.3 cm



"Positioning equipment should be used in a safe manner and according to manufacturer's written instructions."

2009 Perioperative Standards and Recommended Practices. RP: Positioning the Patient, AORN, Inc.





BEST SELLING VALUE PRODUCTS

If your O.R. could have nothing else- consider these basic, versatile products...

Donut Head Pads

- Best seller.
- Variety of sizes for optimum positioning.
- Optional polymer lined dish provides premium protection.



40201	Adult	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40212	Adult Tall	8 OD x 3 ID x 3"	20.5 OD x 7.5 ID x 7.5 cm
40202	Pediatric	5 1/2 OD x 2 1/4 ID x 1 1/4 "	14 OD x 5.7 ID x 3.2 cm
40203	Neo-Natal	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40218	Adult Wide	8 OD x 4 1/2 ID x 1 3/4"	20.5 OD x 11.4 ID x 4.5 cm
40201D	Adult with Center Dish	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40202D	Pediatric with Center Dish	5 1/2 OD x 2 1/4 ID x 1 1/4 "	14 OD x 5.7 ID x 3.2 cm
40203D	Neo-Natal with Center Dish	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm



Action® O.R. Overlay, Small

- Also known as the Hip Pad.
- Most versatile size for maximum utility.
- Use with Model 40100 or 40103 for full table coverage.

40105	20 x 20 x 1/2"	50.8 x 50.8 x 1.3 cm
--------------	----------------	----------------------

Ulnar Nerve Protector



- Protects ulnar nerve while allowing full access to arm.
- Convenient, quick hook-and-loop attachments.

40304	Small	15 3/4 x 6 x 1/2"	40 x 15 x 1.3 cm
40304L	Large	19 x 8 1/2 x 1/2"	48 x 21.6 x 1.3 cm



Roll Covers


- Wraps around staph check-type rolls your facility already has.
- Upgrade basic chest rolls for pressure redistribution.
- Features hook-and-loop attachments.

40600	Chest Roll Cover	18 x 18 x 1/2"	46 x 46 x 1.3 cm
40601	Axillary Roll Cover	17 x 12 x 1/2"	43 x 30.5 x 1.3 cm

Patient Positioners

- Versatile positioners eliminate use of old-fashioned sandbags.
- Polymer positioners are a good solution to unsafe rolled towels, IV bags and sandbags.


40602	Small	10 x 4 1/2 x 2"	25.4 x 11.4 x 5 cm
40605	Large	14 x 4 1/2 x 2"	35.5 x 11.4 x 5 cm
40615	Extra Long	18 x 4 x 2"	46 x 10 x 5 cm



Heel Support

- Secure and protect heel area without hook-and-loop attachment.
- Channel in back of cup cradles patient's Achilles tendon area.

40502	6 x 3 1/4" x 1 3/4"	15 x 8.3 x 4.5 cm
--------------	---------------------	-------------------



Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

KIT002	4 x 4"	10 x 10 cm
KIT003	36 x 4"	91 x 10 cm





Pediatric Positioner Set

- Multipurpose set protects neonates to toddlers during a wide variety of procedures.
- Rectangular pad with indented channel allows chest expansion.
- Pediatric horseshoe pad (Model 40205) supports head in prone or lateral position.
- Two different size chest rolls (Model 40642 and 40646).

40614 Rectangular Pad: 20 x 9 x 2"

50.8 x 22.9 x 5 cm



Candy Cane Stirrup Pad Set

- Minimize risk of skin breakdown due to metal cane contact.
- Protects the peroneal nerve during lithotomy procedures.
- Eliminates tedious wrapping, unwrapping and re-wrapping.
- Use with Action® Heel Protectors, Model 40401.

40501 (Set of 2)

28 x 2" OD

71 x 5 cm OD



Knee Crutch Pads

- Three sizes for maximum versatility.

40402	Small	11 x 9 x 3/8"	27.9 x 22.9 x .95 cm
40403	Medium	17 x 11 x 1/4"	43 x 27.9 x .6 cm
40404	Large	22 x 20 x 1/2"	56 x 50 x 1.3 cm



Armboard Pads

- Protect patient's bicep, elbow, and forearm areas.
- Use on standard armboard table extensions.

40300	Short	15 x 4 x 1/4"	38 x 10 x .6 cm
40301	Long	24 x 4 1/2 x 3/8"	61 x 11.5 x .95 cm
40309	Wide Short	15 x 6 1/2 x 1/4"	38 x 16.5 x .6 cm
40308	Wide Long	24 x 6 1/2 x 3/8"	61 x 16.5 x .95 cm



Warranty

Action® pads are warranted to be free of defective materials and workmanship at the time when first purchased and are serviceable for a period of two years provided the product is not subject to misuse, negligence, accident or abuse. Accessories are warranted to be free of defects in workmanship at the time of initial purchase only. Specific warranty questions should be directed to Customer Service at 301.797.1414 (Monday - Friday 8:30AM - 5:00PM EST) or your authorized dealer/distributor.

CUSTOM PADS AND POSITIONERS ARE OUR SPECIALTY.

Please call or fax a diagram with your specific requirements. Custom orders cannot be returned. All pads and positioners are imaging compatible and easily disinfected.

TERMS:

Net 30 days, 2% discount on merchandise only for payment within 10 days of invoice date. FOB Hagerstown, MD. Service Charge of 1.5% per month on unpaid balance. Prices may change without notice. Handling charge of \$5.00 per order. \$3.00 drop ship fee may apply.

Action Products, Inc.

954 Sweeney Drive
Hagerstown, MD 21740, USA
800.228.7763
301.797.1414
301.733.2073 (Fax)
service@actionproducts.com
www.actionproducts.com

