





- 4-way stretch fabric for shear reduction
- High-density foam base
- Radiolucent
- Latex Free
- Fluid proof, welded seams
- 100% no-stitch hook and loop attachment
- Easy to clean
- Great pad for heavy frame use
- Nurse friendly ergonomic design
- Other models available upon request

NEW! Action® Original FT Replacement Table Pads are a terrific value. Simple pressure management is maintained with a 4-way stretch fabric atop a 3/8" Akton® polymer and high density foam. The Action® standard infection control design and bacteria resistant fabric create a durable and reliable table pad suitable for any budget.



3" Original FT Pad
57" Male Subject, 170 lbs. Average Pressure (mmHg): 20.58 Room Temp: 70 degrees Fahrenheit

Action® Original FT Table Pads

Models	2 3/8" Original FT- 3-pc. Set	3 3/8" Original FT- 3-pc. Set	2 3/8" Original FT Silver 3-pc. Set	3 3/8" Original FT Silver 3-pc. Set
Steris Amsco 3080/3085	75120	77120	75120DS	77120DS
Steris C-Max	75130	77130	75130DS	77130DS
Skytron 3500	75310	77310	75310DS	77310DS
Skytron 5001, 6001, 6002, 6500 & 6600	75320	77320	75320DS	77320DS
Skytron 6700, 6701	75330	77330	75330DS	77330DS
Skytron 3501, 3600, 3600B	75370	77370	75370DS	77370DS
Flat Armboard Pads	79100	79600	79100DS	79600DS

Visit us at www.4yourOR.com and join the pressure management discussion on our social network:
<http://blog.actionproducts.com>, and .

