




- Shear Smart® 360° stretch surface
- High-density foam base
- 4-way stretch fabric for shear reduction
- Radiolucent
- Latex Free
- Fluid proof, welded seams
- 100% no-stitch hook and loop attachment
- Easy to clean
- Nurse friendly ergonomic design
- Other models available upon request


The Original pads are like two products in one: a simple, premium construction table pad and an Action® Overlay surface. Action invented the Shear Smart® polymer surface to protect patients from O.R. acquired pressure ulcers and has taken its flagship pad to a higher level. The integrated 3/8" Akton® polymer overlay evenly redistributes heat and weight, prevents bottoming out, minimizes slippage, and improves the longevity of a standard foam table pad. The Original is now available in both 2 3/8" and 3 3/8" thicknesses. This unique combination satisfies the needs of the patient, the O.R. staff and the budget-driven facility.



V Value, Versatility, Clinical Effectiveness--- All this from the pad you have trusted for over 15 years for the best shear and pressure management benefits in the industry. Who else can say that?

Action® Original Table Pads

Models	2 3/8" Original 3-pc. Set	3 3/8" Original 3-pc. Set
Steris Amsco 3080/3085	70120	72120
Steris C-Max	70130	72130
Skytron 3500	70310	72310
Skytron 5001, 6001, 6002, 6500 & 6600	70320	72320
Skytron 6700, 6701	70330	72330
Skytron 3501	70370	72370
Flat Armboard Pads	74100	74600



3" Original Pad
60" Male Subject, 180 lbs. Average Pressure (mmHg): 19.4 Room Temp: 70 degrees Fahrenheit

Visit us at www.4yourOR.com and join the pressure management discussion on our social network:

<http://blog.actionproducts.com>, and .

