



- Clinically proven
- Redistributes heat and weight
- Soft shear-resistant film
- Will not support bacterial growth
- Biocompatible
- Fast and easy cleaning
- Latex free
- Does not leak or flow
- Tissue equivalent

Introducing the Action® Orthopedic Hip Positioner...

The Action® Orthopedic Hip Positioner was originally commissioned by Dr. William Hozack of the renowned *Rothman Institute Orthopaedics* in Philadelphia, PA for use during total hip replacements. Designed to complement the Anterior Supine Approach to THR, this unique positioner, often called a “bump”, was engineered to provide consistent, optimal site access and patient stability while protecting the patient’s skin integrity. Our soft Akton® polymer encases a rigid foam core to create a critical frame-like stability while simultaneously redistributing pressure and shear forces that can damage the sacral tissue.

- Improves patient safety by redistributing pressure and minimizing shear
- Eliminate use of potentially harmful makeshift items such as rolled/ taped blankets, towels, IV bags, etc...
- Weighted platform provides rigid, secure base required during reaming
- Reusable and durable, easy to clean
- Physician preferred item

	Dimensions		Weight	
41740	16 x 8.5 x 2.63"	40.64 x 21.59 x 6.68 cm	5.56 lbs.	2.52 kg



“We have used it exclusively for our hip replacement surgeries for several years without complication and today, we even train our fellows to perform the surgeries using this unique positioner.”

Dr. William Hozack,
Rothman Institute Orthopaedics

Introducing the Peg Board Overlay and Post Pad Set...

Protect your patients from damaging pressure and shear forces with our newest Action® pad set designed to cushion the standard peg board device. The Action® Overlay and Post Pad Set protect skin integrity in the high-pressure lateral position, by redistributing pressure on fragile bony prominences. The post pads have a hook and loop attachment to universally fit your morph or peg board posts.

The pads are sold as a set or replacement pads may be purchased individually. Peg board and posts are not included.



SKU# 40719 Set includes:

		Dimensions		Weight	
40719 A	4 Short Post Pad	5 x 5.5 x 0.5"	12.7 x 14 x 1.3 cm	0.4 lbs.	0.2 kg
40719 B	4 Tall Post Pad	5 x 8.5 x 0.5"	12.7 x 21.6 x 1.3 cm	0.6 lbs.	0.3 kg
40100	Action® Overlay Pad	46 x 20 x 0.5"	116.8 x 50.8 x 1.3 cm	19 lbs.	8.6 kg

