


- 4-way stretch antimicrobial fabric surface.
- Continuous layer of 3/8" Akton® polymer.
- Soft viscofoam mid-layer.
- High-density foam base.
- Fluid proof, welded seams.
- 100% no-stitch hook-and-loop attachment.
- Proprietary hidden vent system.
- Easy to clean waterfall construction.
- Available in multiple thicknesses.
- Armboards available.
- Custom pads are available. Call for details.



3" Performance Pad
 5'11" Male Subject, 175 lbs. Average Pressure (mmHg): 18.72 Room Temp: 68 degrees Fahrenheit
 Equipment: FSA Pressure Map System Range: 1 – 150 mmHg

Introducing our **Performance** and **Performance PLUS Series** of table pads that will help you achieve your zero tolerance goals for HAPUs (Hospital Acquired Pressure Ulcers) and fall well within your budget.


Our pads are clinically engineered to give the maximum in decubiti prevention, infection control, and stability for the ultimate in patient safety. With two styles and multiple thicknesses, we have the perfect mattress for any budgetary requirement.


Action® Performance Series O.R. Table Pads

55120	2" Fabric Top	Fits Amsco 3080/3085
55130	2" Fabric Top	Fits Steris C-Max
55310	2" Fabric Top	Fits Skytron 3500
55320	2" Fabric Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
55330	2" Fabric Top	Fits Skytron 6700/ 6701
55370	2" Fabric Top	Fits Skytron 3501
57120	3" Fabric Top	Fits Amsco 3080/3085
57130	3" Fabric Top	Fits Steris C-Max
57310	3" Fabric Top	Fits Skytron 3500
57320	3" Fabric Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
57330	3" Fabric Top	Fits Skytron 6700/ 6701
57370	3" Fabric Top	Fits Skytron 3501
59100	2" Fabric Top	Universal Armboard
59600	3" Fabric Top	Universal Armboard




Durable, welded, midline seam and proprietary hidden vent structure enhance infection control for maximum pad life.





- Shear Smart® surface of exposed 3/8" Akton® polymer.
- Multi-directional stretch film covering on surface.
- Soft viscofoam mid-layer.
- High-density foam base.
- Tissue equivalent.
- Does not support bacterial growth.
- Fluid proof, welded seams.
- 100% no-stitch hook-and-loop attachment.
- Proprietary hidden vent system.
- Easy to clean waterfall construction.
- Available in multiple thicknesses.
- Armboards available.
- Custom pads are available. Call for details.



3" Performance Plus Pad

5' 11" Male Subject, 175 lbs. Average Pressure (mmHg): 22.33 Room Temp: 69 degrees Fahrenheit

Equipment: FSA Pressure Map System Range: 1 – 150 mmHg

Action® Performance Plus Series O.R. Table Pads

50120	2" Polymer Top	Fits Amsco 3080/3085
50130	2" Polymer Top	Fits Steris C-Max
50310	2" Polymer Top	Fits Skytron 3500
50320	2" Polymer Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
50330	2" Polymer Top	Fits Skytron 6700/ 6701
50370	2" Polymer Top	Fits Skytron 3501
52120	3" Polymer Top	Fits Amsco 3080/3085
52130	3" Polymer Top	Fits Steris C-Max
52310	3" Polymer Top	Fits Skytron 3500
52320	3" Polymer Top	Fits Skytron 5001/ 6001/ 6002B/ 6500/ 6600
52330	3" Polymer Top	Fits Skytron 6700/ 6701
52370	3" Polymer Top	Fits Skytron 3501
54100	2" Polymer Top	Universal Armboard
54600	3" Polymer Top	Universal Armboard



AKTON polymer in a proprietary soft film that reduces shear effects for optimal patient safety.



Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

KIT002 4 x 4" 10 x 10 cm

KIT003 36 x 4" 91 x 10 cm

- VALUE product essential for any operating room.

