

- Cited by AORN best practices for Pressure Ulcer Prevention.
- Provide pressure redistribution and reduce shear effects across entire table surface.
- Minimize patient movement.
- Low profile, simple design.
- ½" Akton® polymer for maximum effectiveness.



## Hourglass Hand Table Overlay

- Appropriate for lengthy micro-hand procedures.

**40330** 32 x 24 x ½" 81.3 x 61 x 1.2 cm



## Mayo Stand Cover

- Use for orthopedic procedures.
- Cushion costly instruments.
- Pillowcase style with rubber bottom.

**40717** 21 x 18 ½ x ¾" 53.3 x 47 x 1.6 cm

Pillowcase



Fitted



## Action® Disposable Overlay Covers

- Infection control product.
- Fitted sheet and pillowcase styles.
- Extend the service life of the pad.
- Improves room turnover rate.
- Case quantities available.

### Pillowcase Style

**40102** Made for table pad models 40100 & 40101 Case of 25

### Fitted Sheet Style

**40800** Made for table pad model 40100 Case of 25

**40801** Made for table pad model 40101 Case of 25

**40803** Made for table pad model 40103 Case of 25

**40804** Made for table pad model 40104 Set of 3

**40805** Made for table pad model 40105 Each



## Action® O.R. Overlay with Perineal Cutout

- Ideal for lithotomy position.
- Wide range of applications in both regular and cysto operating rooms.

## Action® O.R. Overlay, 3-Piece Set

- For use with segmented tables.
- Includes head, torso, and foot section pads.
- Provides pressure and shear reduction across entire table surface.

## Action® O.R. Overlays

- Standard sizes and customs available.

<b>40101</b>	Large (Full-length)	72 x 20 x ½"	183 x 50 x 1.3 cm
<b>40100</b>	Medium	46 x 20 x ½"	117 x 50 x 1.3 cm
<b>40105</b>	Small - Hip pad	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm <input type="checkbox"/>
<b>40103</b>	Medium with cutout	46 x 20 x ½" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
<b>40104A</b>	Head	9 x 20 x ½"	22.9 x 50.8 x 1.3 cm
<b>40104B</b>	Torso	39 ½ x 20 x ½"	100.3 x 50.8 x 1.3 cm w/ perineal cutout
<b>40104C</b>	Foot	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm w/ perineal cutout



## Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any polymer surface.

**KIT002** 4 x 4" 10 x 10 cm **KIT003** 36 x 4" 91 x 10 cm

- VALUE product essential for any operating room.

