



Pediatric Chest Rolls

40639	6 x 2"	15 x 5 cm
40640	8 x 2"	20 x 5 cm
40641	10 x 2"	25.5 x 5 cm
40642	12 x 2"	30.5 x 5 cm
40643	6 x 3"	15 x 7.5 cm
40644	8 x 3"	20 x 7.5 cm
40645	10 x 3"	25.5 x 7.5 cm
40646	12 x 3"	30.5 x 7.5 cm



Pediatric / Neo-Natal Positioners

40648	6 x 1½ x ¾"	15 x 3.8 x 2 cm
40649	6 x 1 x 1"	15 x 2.5 x 2.5 cm
40650	6 x 1½ x 1½"	15 x 3.8 x 3.8 cm



Pediatric Positioner Set V

- Multi-purpose set protects neonates to toddlers during a wide variety of procedures
- Rectangular pad with indented channel allows chest expansion
- Pediatric Horseshoe Pad (40205) supports head in prone or lateral position
- Includes two different size Chest Rolls (40642 and 40646)

40614	Four piece pediatric positioning set
40647	Pad w/ channel: 20 x 9 x 2" 50.8 x 22.9 x 5 cm



40202



40205



40221



40217P

Pediatric Headrests

40202	Pediatric Donut Head Pad	5½ OD x 2¼ ID x 1¼"	14 OD x 5.7 ID x 3.2 cm
40219	Pediatric/ Neo-Natal Donut Head Pad	4 OD x 2 ID x 1"	10.2 OD x 6.4 ID x 2.5 cm
40203	Neo-Natal Donut Head Pad	3¼ OD x 1½ ID x ¾"	8.3 OD x 3.8 ID x 2 cm
40202D	Pediatric Donut Head Pad with Center Dish	5½ OD x 2¼ ID x 1¼"	14 OD x 5.7 ID x 3.2 cm
40203D	Neo-Natal Donut Head Pad with Center Dish	3¼ OD x 1½ ID x ¾"	8.3 OD x 3.8 ID x 2 cm
40205	Pediatric Horseshoe Head Pad	5½ OD x 2¼ ID x 1¼"	14 OD x 5.7 ID x 3.2 cm
40220	Pediatric/ Neo-Natal Horseshoe Head Pad	4 OD x 2 ID x 1"	10.2 OD x 6.4 ID x 2.5 cm
40206	Neo-Natal Horseshoe Head Pad	3¼ OD x 1½ ID x ¾"	8.3 OD x 3.8 ID x 2 cm
40217P	Ophthalmic Cradle Headrests Pediatric	7 ½ x 8 ½ x 2 ¼"	19 x 21.6 x 5.7 cm
NEW 40221	Pediatric Head Pad with Center Dish	6 x 6 x 1"	15.2 x 15.2 x 2.5 cm

V - VALUE product essential for any operating room.





Forehead Protector

- Protects infant, child, or adult patient from pressure exerted by oxygen tubes, etc.
- Semi-circular with flat bottom

40208	2½ x 1¾ x ¾"	6.3 x 4 x 1.6 cm
--------------	--------------	------------------



Isolette Pad

- Consistent minimization of pressure and shear in NIC units to O.R.

40120L	Isolette Pad	26 x 16 x ½"	66 x 40.6 x 1.3 cm
---------------	--------------	--------------	--------------------

Also available:

40120	Bassinet Pad	26 x 12 x ½"	66 x 30.5 x 1.3 cm
--------------	--------------	--------------	--------------------



Flat



Rolled

Pediatric Ulnar Nerve Protector

- Protects ulnar nerve while allowing full access to arm
- Convenient, quick hook-and-loop attachments

NEW	40304P	Pediatric	10 x 3 x ¼"	25.4 x 7.6 x .6 cm
	40304	Small	15¾ x 6 x ½"	40 x 15 x 1.3 cm




Patient Positioners

- Versatile positioners eliminate use of old-fashioned sandbags.
- Polymer positioners are a good solution to unsafe rolled towels, IV bags and sandbags.

40602	Small	10 x 4 ½ x 2"	25.4 x 11.4 x 5 cm
40605	Large	14 x 4 ½ x 2"	35.5 x 11.4 x 5 cm
40615	Extra Long	18 x 4 x 2"	46 x 10 x 5 cm

Universal Spot Protectors



- Adaptive use padding options for those hard-to-protect areas unique to each patient
- Designed to enhance pressure redistribution and improve shear performance

40405	Mini	1½ OD x ¼"	3.8 OD x .6 cm	Pkg of 12
40406	Small	3¼ OD x ¼"	8.3 OD x .6 cm	Pkg of 12
40407	Large	5¼ OD x ¼"	13.3 OD x .6 cm	Pkg of 6

- VALUE product essential for any operating room.

