


ACTION® Disposable Covers

- ✓ Soft 360° stretch film
- ✓ Easy on and off design
- ✓ Helps reduce clean up time
- ✓ Extends the service life



Pillowcase style--
full length overlay cover

Disposable Overlay Covers

Action's waterproof, non-sterile Disposable Overlay Covers are made of a soft, flexible film and is available in two styles. Designed to help minimize exposure to fluids and to help extend the life of the overlay. The "pillowcase" style fully encases the overlay. The "fitted sheet" style features pockets at both ends of the pad and has a 4" film drape along the length of the overlay.

Fitted Sheet Style		
40800	Made for overlay model 40100- torso overlay	Case of 25
40801	Made for overlay model 40101- full length pad	Case of 25
40803	Made for overlay model 40103- torso with cut out	Case of 25
40804	Made for overlay model 40104	3 Piece Set
40805	Made for small hip pad	Each
Pillowcase Style		
40102	Made for table pad models 40100 & 40101	Case of 25

Fitted Sheet style--torso with cut out



Disposable Chest Roll Covers

Action's waterproof, non-sterile Disposable Chest Roll Covers are made of a soft, flexible film designed to extend the life of the chest roll and minimize exposure to fluids during use. The slip-on covers mimic size and shape of the chest roll and have an overlapping film closure along the flat bottom.

Purchase in cases of 25.

40616DC	12 x 4"	30.5 x 10 cm
40617DC	16 x 4"	40.6 x 10 cm
40618DC	20 x 4"	50.8 x 10 cm
40619DC	12 x 6"	30.5 x 15 cm
40620DC	16 x 6"	40.6 x 15 cm
40621DC	20 x 6"	50.8 x 15 cm



Action Products, Inc.
 800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2021 Action Products, Inc.



CERTIFIED