



## OPERATION DryMax<sup>®</sup>

- Single-use disposable product
- Safer work environment reduces slips, trips and falls
- Improves clean up efficiency for faster turnover times
- Improves infection control by trapping and retaining contaminated fluids
- Reduce red bag costs
- Lower laundry bills from cleaning towels and blankets used to soak up liquids
- Perfect fluid management for operating room, procedure or sink areas
- Great for Orthopedics, Urology, Labor & Delivery, Central Sterilization and more
- **Adheres to OSHA, JCAHO and CMS standard regulations**




### DryMax 2.4 - superabsorbent

- Traps and retains fluids
- Low profile, ultra-thin mat
- Absorbs from top and bottom
- Safer work environment
- Fast, efficient turnover
- **9 liter capacity**




### DryMax TRIPLE - superabsorbent

- Traps and retains fluids
- Flexible 3-part mat
- Full OR table length
- Low profile, ultra-thin mat
- Absorbs from top and bottom
- **15 liter capacity**



### DryMax XL - absorbent

- Anti-skid design improves safety
- Large surface area keeps floors dry
- Can be placed under a stool
- Ideal to use under sink area
- Good for endoscopy based surgery and urology cases



### DryMax COMBIMAT - flexible absorbent

- High absorbency combination of the 2.4 and XL
- Low profile, ultra-thin mat
- Anti-skid barrier
- Green mat absorbs from top and bottom
- Can easily be folded to protect shoes and floors
- Good for endoscopy based surgery, gynecology and urology cases

|                 | Model   | Size         | Weight | Unique Quality              | Packaging            | Absorbency:  = 1 liter            |
|-----------------|---------|--------------|--------|-----------------------------|----------------------|-----------------------------------|
| DryMax 2.4      | DM20241 | 28 x 14 1/2" | 120 g  | Absorbs from top and bottom | 25 pcs/bag, 100/case | 9 lt of water<br>3.5 lt of NaCl   |
| DryMax TRIPLE   | DM20243 | 85 x 14 1/2" | 200 g  | Folds around objects        | 40/case              | 15 lt of water<br>6 lt of NaCl    |
| DryMax XL       | DM22201 | 39 x 29 1/2" | 160 g  | Non-skid                    | 10 pcs/bag, 50/case  | 1.5 lt of water<br>1.5 lt of NaCl |
| DryMax COMBIMAT | DM22202 | 39 x 43"     | 350 g  | 2 mats in 1 XL + 2.4        | 5 pcs/bag, 30/case   | 15 lt of water<br>7 lt of NaCl    |