

ACTION® OVERLAYS



Overlays

The Action® O.R. Overlay is a versatile polymer pad used for a wide range of procedures. Its shear reducing qualities combined with the ability to distribute heat and pressure make it a valuable pressure management tool. Used in longer procedures or for patients with fragile skin integrity, the overlay is an affordable upgrade for basic foam table pads that provides enhanced patient protection. The overlays are portable and can be exchanged between rooms as needed and are available in all standard OR table footprints in 1/2" thickness. Disposable covers are available for select overlay models.


40105	Small	20 x 20 x 1/2"	50.8 x 50.8 x 1.3 cm
40100	Medium	46 x 20 x 1/2"	117 x 50 x 1.3 cm
40109	Medium Large	63 x 20 x 1/2"	160 x 50 x 1.3 cm
40101	Full Length	72 x 20 x 1/2"	183 x 50 x 1.3 cm
40103	Perineal Cutout	46 x 20 x 1/2" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
40104	3-piece Set	See A, B and C	
40104A	Head	9 x 20 x 1/2"	22.9 x 50.8 x 1.3 cm
40104B	Torso w/ perineal cutout	39 1/2 x 20 x 1/2"	100.3 x 50.8 x 1.3 cm
40104C	Foot w/ perineal cutout	20 x 20 x 1/2"	50.8 x 50.8 x 1.3 cm



Universal Spot Protectors

This series of thin, round pads is designed to provide adaptive use padding options for those hard-to-protect areas unique to each patient.

40405	Mini	1 5/8" D x 1/4"	4.1 D x .6 cm	Pkg of 12
40406	Small	3 1/4" D x 1/4"	8.3 D x .6 cm	Pkg of 12
40407	Large	5 1/2" D x 1/4"	14 D x .6 cm	Pkg of 6



AORN 2023 Guidelines, Pressure Injury Prevention

Selecting and Using Support Surfaces, section 3.2

“Position patients on surfaces that reduce the potential for PI by redistributing pressure and reducing the risk of unintended patient movement.... Static surfaces such as viscoelastic polymer (i.e., gel) overlays are designed to reduce shearing and to support the patient's weight without becoming fully compressed under the weight of the patient's body.^{5,46,74} Gel overlays are radiolucent, latex-free, and reusable.⁷⁴

Advantages of AKTON® Polymer



- Clinically proven to help prevent pressure injuries
- Distributes weight evenly
- Distributes heat evenly
- Will not leak or flow if punctured
- Will not bottom out
- Will not support bacterial growth
- Imaging compatible
- Biocompatible
- Reusable
- Cost effective
- Tissue equivalent with 360 degree stretch film to help reduce shear
- Hypo-allergenic
- Not made with natural rubber latex
- Free of phthalates, plasticizers, and silicone

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.

ACTION products
feature our proprietary
Akton® polymer





Mayo Stand Cover

Use to prevent pressure and nerve damage when limb is placed on the Mayo Stand for positioning during orthopedic procedures.

40717 21 x 18 1/2 x 5/8" 53.3 x 47 x 1.6 cm



Hourglass Hand Table Overlay

Designed to be used in conjunction with an Hour Glass Hand Table, this product is particularly appropriate for lengthy micro-hand procedures that require additional pressure relieving protection.

40330 32 x 24 x 16 x 1/2" 81.3 x 61 x 40.6 x 1.2 cm

Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any Akton® polymer surface.


See *Care and Use Guide (IFU)* at www.actionproducts.com for information on resurfacing our Akton® polymer.

KIT002 4 x 4" 10 x 10 cm

KIT003 36 x 4" 91.4 x 10 cm

Disposable Overlay Covers

Action's waterproof, non-sterile Disposable Overlay Covers are made of a soft, flexible film and are available in two styles. Designed to help minimize exposure to fluids and to help extend the life of the overlay. The "pillowcase" style fully encases the overlay. The "fitted sheet" style features pockets at both ends of the pad and has a 4" film drape along the length of the overlay.



Pillowcase style--
full length overlay cover

Fitted Sheet Style

40800	Made for overlay model 40100- torso overlay	Case of 25
40801	Made for overlay model 40101- full length pad	Case of 25
40803	Made for overlay model 40103- torso with cut out	Case of 25
40805	Made for small hip pad	Case of 25

Pillowcase Style

40102	Fits overlay models 40100 & 40101	Case of 25
--------------	-----------------------------------	------------



Fitted Sheet style--torso with cut out

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073

www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.



M3001-001