


## Overlays

The Action® O.R. Overlay is a versatile polymer pad used for a wide range of procedures. Its shear reducing qualities combined with the ability to distribute heat and pressure make it a valuable pressure management tool.

Used in longer procedures or for patients with fragile skin integrity, the overlay is an affordable upgrade for basic foam table pads that provides enhanced patient protection. The overlays are portable and can be exchanged between rooms as needed and are available in all standard OR table footprints in 1/2" thickness. Disposable covers are available for select overlay models.

<b>40105</b>	Small	20 L x 20 W x ½ H"	50.8 x 50.8 x 1.3 cm
<b>40100</b>	Medium	46 L x 20 W x ½ H"	117 x 50 x 1.3 cm
<b>40109</b>	Medium Large	63 L x 20 W x ½ H"	160 x 50 x 1.3 cm
<b>40101</b>	Full Length	72 L x 20 W x ½ H"	183 x 50 x 1.3 cm
<b>40103</b>	Perineal Cutout	46 L x 20 W x ½ H" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
<b>40104</b>	3-piece Set	See A, B and C	
<b>40104A</b>	Head	9 L x 20 W x ½ H"	22.9 x 50.8 x 1.3 cm
<b>40104B</b>	Torso w/ perineal cutout	39 ½ L x 20 W x ½ H"	100.3 x 50.8 x 1.3 cm
<b>40104C</b>	Foot w/ perineal cutout	20 L x 20 W x ½ H"	50.8 x 50.8 x 1.3 cm



## AORN 2024 Guidelines, Pressure Injury Prevention Selecting and Using Support Surfaces, section 3.2

“Position patients on surfaces that reduce the potential for PI by redistributing pressure and reducing the risk of unintended patient movement.... Static surfaces such as viscoelastic polymer (i.e., gel) overlays are designed to reduce shearing and to support the patient's weight without becoming fully compressed under the weight of the patient's body.<sup>5,46,74</sup> Gel overlays are radiolucent, latex-free, and reusable.<sup>74</sup>


## Disposable Overlay Covers

Action's fluid-resistant, non-sterile Disposable Overlay Covers are made of a soft, flexible film and are available in two styles. They are designed to help reduce exposure to fluids and to help extend the life of the overlay.

The “pillowcase” style is closed at one end and fully encases the overlay like a sheath.

The “fitted sheet” style covers the top of the pad and is held in place by two pockets at each end of the cover. This cover has a 4" film drape that runs the length of the cover on both sides.

Fitted Sheet Style		
<b>40800</b>	Made for overlay model 40100- medium overlay	Case of 25
<b>40801</b>	Made for overlay model 40101- full length overlay	Case of 25
<b>40803</b>	Made for overlay model 40103- medium with cut out	Case of 25
<b>40805</b>	Made for overlay model 40105- small overlay	Case of 25
Pillowcase Style		
<b>40102</b>	Fits overlay models 40100 & 40101	Case of 25



## Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073  
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products  
feature our proprietary  
Akton® polymer





### Mayo Stand Cover

Use to prevent pressure and nerve damage when limb is positioned on the Mayo Stand during procedures. The 5/8" polymer pad is sealed along three sides to a thick, durable film creating a pocket that slides easily on to the top of the Mayo stand.


**40717** 21 L x 18 1/2 W x 5/8 H" 53.3 x 47 x 1.6 cm



### Hourglass Hand Table Overlay


This reusable 1/2" polymer overlay designed for surgical procedures involving the forearm, elbow, or hand utilizing the hourglass hand table. Good for lengthy hand microsurgeries requiring additional pressure and shear relieving support. Fits most standard hourglass hand tables.

**40330** 32 L x 24 W x 1/2 H" 81.28 x 60.96 x 1.27 cm  
16"W middle section 40.6 cm middle section



### Universal Spot Protectors

This series of thin, round pads is designed to provide adaptive use padding options for those hard-to-protect areas unique to each patient. The shape is designed to enhance pressure redistribution and improve shear performance.



<b>40405</b>	Mini	1 5/8 D x 1/4 H"	4.1 D x .6 cm	Pkg of 12
<b>40406</b>	Small	3 1/4 D x 1/4 H"	8.3 D x .6 cm	Pkg of 12
<b>40407</b>	Large	5 1/2 D x 1/4 H"	14 D x .6 cm	Pkg of 6

### NICU Pads

Provides pressure and shear management surface for pediatric and neo-natal patients.

<b>40120L</b>	Isolette Pad	26 L x 16 W x 1/2 H"	66 x 40.6 x 1.3 cm
<b>40121</b>	Isolette Pad	29 L x 19 W x 1/2 H"	74 x 48.2 x 1.3 cm
<b>40120</b>	Bassinet Pad	26 L x 12 W x 1/2 H"	66 x 30.5 x 1.3 cm




### Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any Akton® polymer surface.

See Care and Use Guide (IFU) at [www.actionproducts.com](http://www.actionproducts.com) for information on resurfacing our Akton® polymer.

<b>KIT002</b>	4 L x 4 W"	10 x 10 cm
<b>KIT003</b>	36 L x 4 W"	91.4 x 10 cm



### Advantages of AKTON® Polymer

- Excellent shear reducing qualities
- Soft, flexible film
- Enhances patient stability
- Evenly distributes heat and pressure
- Clinically proven to help prevent pressure injuries
- Reusable & Cost Effective
- Does not support bacterial growth
- Hypoallergenic
- Imaging compatible
- Biocompatible
- Hypo-allergenic
- Not made with natural rubber latex
- Free of phthalates, plasticizers & silicone
- Free of DEHP

