




# CHEST / TORSO / HIP POSITIONING PRODUCTS



Adult		
40616	12 ½ L x 4 ¾ W x 3 ¾" H	31.75 x 11.11 x 8.57 cm
40617	16 ½ L x 4 ¾ W x 3 ¾" H	41.91 x 11.11 x 8.57 cm
40618	20 ½ L x 4 ¾ W x 3 ¾" H	52.07 x 11.11 x 8.57 cm
40619	12 L x 6 W x 5" H	30.48 x 15.24 x 12.7 cm
40620	16 L x 6 W x 5" H	40.64 x 15.24 x 12.7 cm
40621	20 L x 6 W x 5" H	50.8 x 15.24 x 12.7 cm

## Chest Rolls

This line of flat bottom chest rolls is designed to provide stable support in the prone position. They are available in a wide range of heights and lengths to accommodate virtually any patient, adult to neonate.



Pediatric		
40639	6 L x 2 W x 2" H	15.24 x 5.08 x 5.08 cm
40640	8 L x 2 W x 2" H	20.32 x 5.08 x 5.08 cm
40641	10 L x 2 W x 2" H	25.4 x 5.08 x 5.08 cm
40642	12 L x 2 W x 2" H	30.48 x 5.08 x 5.08 cm
40643	6 L x 3 W x 3" H	15.24 x 7.62 x 7.62cm
40644	8 L x 3 W x 3" H	20.32 x 7.62 x 7.62cm
40645	10 L x 3 W x 3" H	25.4 x 7.62 x 7.62cm
40646	12 L x 3 W x 3" H	30.48 x 7.62 x 7.62cm

## Disposable Chest Roll Covers


Action's fluid-resistant, non-sterile Disposable Chest Roll Covers are made of a soft, flexible film designed to extend the life of the chest roll and reduce exposure to fluids during use.

The slip-on covers mimic size and shape of the chest roll and have an overlapping film flap along the flat bottom.

*Purchase in cases of 25.*



Fits Action® Chest Roll	
40616DC	40616
40617DC	40617
40618DC	40618
40619DC	40619
40620DC	40620
40621DC	40621



## Chest Positioning Set

These positioners have a contoured design patterned to uniformly conform to a patient's shoulder to chest area, with shaped edges for the axilla. They provide upper body positioning while alleviating pressure on the thorax area.

40625	Left	28 - 23 ¾ L x 8 x 3 x 2 ½" H	71.12 x 58.42 x 20.32 x 7.62 x 6.35 cm
40626	Right	28 - 23 ¾ L x 8 x 3 x 2 ½" H	71.12 x 58.42 x 20.32 x 7.62 x 6.35 cm

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073  
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products  
feature our proprietary  
Akton® polymer





# ACTION® CHEST / TORSO / HIP POSITIONING PRODUCTS

## Lateral Positioner

Designed to protect the patient while in the lateral position. Accommodates both shoulder and upper arm. Provides pressure relief and nerve protection. Includes a lightweight foam base with AKTON polymer pad overlay.


**40610** 28 L x 19 1/2 W x 6 - 1" H 71.12 x 49.53 x 15.24 - 2.54 cm



## Patient Positioner

AORN encourages the use of Axillary Rolls in place of sandbags or IV bags. These extremely versatile positioners are good for many procedures.

<b>40602</b> Small	10 L x 4 1/4 W x 1 7/8" H	25.4 x 10.8 x 4.76 cm
<b>40605</b> Long	14 L x 4 W x 2" H	35.56 x 10.16 x 5.08 cm
<b>40615</b> Extra Long	18 L x 4 W x 2" H	45.72 x 10.16 x 5.08 cm



## Graduated Thoracic Positioner


The graduated, or sloping, configuration of this pad encourages full chest extension. It is designed for cardiovascular or thoracic procedures, which mandate full chest hyperextension.

**40604** 14 L x 9 W x 2 - 1 1/2" H 35.56 x 22.86 x 5.08 - 3.81 cm

## Axillary Roll Cover

Designed to upgrade basic foam positioning rolls with reusable polymer wrap-around covers. Covers are adjustable and are secured by hook and loop attachment.

**40601** 17 L x 12 W x 1/2" H 43.18 x 30.48 x 1.27 cm



## Peg Board Pad Set

The Action® Peg Board Pad set complements most standard peg board positioning systems. This multi-piece set is designed to protect patient skin and soft tissue from damaging pressure and shear forces while in the high-pressure lateral position.

The medium-length overlay (Item #40100) lies atop the peg board (beneath the patient) while the adjustable peg post pads of varying lengths easily wrap around the pegs and are secured by hook and loop attachment. Overlay and peg pads are available as a set, or as individual units.

Set (#40719) includes:

- 1 Medium Overlay Pad, #40100
- 4 Short Peg Post Pads, #40719A
- 4 Tall Peg Post Pads, #40719B

*Peg board and peg posts not included*

<b>40719</b>	Peg Board and Post Pad Set	Set includes (1) 40100, (4) 40719A, (4) 40719B
<b>40719A</b>	Short Post Pad	5 L x 5 W x 5 1/2" H 12.7 x 12.7 x 13.97 cm
<b>40719B</b>	Tall Post Pad	5 L x 5 W x 8 1/2" H 12.7 x 12.7 x 20.32 cm
<b>40100</b>	Medium Overlay	46 L x 20 W x 1/2" H 117 x 50 x 1.3 cm

## Chest Roll Cover

Designed to upgrade basic foam chest rolls with reusable polymer wrap-around covers. Covers are adjustable and are secured by hook and loop attachment.

**40600** 18 L x 18 W x 1/2" H 45.72 x 45.72 x 1.27 cm



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073  
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.



M3006-002



# CHEST / TORSO / HIP POSITIONING PRODUCTS



## Orthopedic Hip Positioner

Designed to complement the Anterior Supine Approach to THR, this unique hip positioner was engineered to provide consistent, optimal site access and patient stability while protecting the patients skin integrity.

**41740** 8 1/2 L x 16 W x 2.63" H 21.59 x 40.64 x 6.68 cm



## Shoulder Brace, 2-Piece Set

Reusable 1" thick polymer pads designed for basic shoulder braces. Adjustable pads wrap around braces and are secured by hook and loop attachments. Sold as a pair.

**40302** 5 L x 6 W x 1" H 12.7 x 15.24 x 2.54 cm

## Action® Comfort Stool Cover

Action® Comfort Stool Cover - 3/8" thick Akton® polymer pad is bonded to a 4" nylon oxford band with drawstring closure.



**40900** 13 D x 3/8" H 33.02 x .95 cm

**40901** 14 D x 3/8" H 35.56 x .95 cm

**40902** 15 D x 3/8" H 38.1 x .95 cm

## Resurfacing Kits

Quick and easy resurfacing of minor cuts or abrasions on any Akton® polymer surface.



**KIT002** 4 x 4" 10.16 x 10.16 cm

**KIT003** 36 x 4" 91.44 x 10.16 cm



## Positioning Strap

Welded seams and fluid-resistant film combined with a soft patient-facing polymer surface make this a must have for the operating room. Radiolucent and MRI safe.

Hook and loop attachments secure the strap to standard tables. The longer strap was specifically designed for larger patients. *This strap is not meant to restrain patient.*



<b>40713</b>	57 L x 4 W x 1/4" H	144.78 x 10.16 x .64 cm	Adjustable range: 28 x 43"	71.12 x 109.22 cm
	AKTON polymer surface measures 98 x 8.2 x 0.6cm			
<b>40713L</b>	90 L x 4 W x 1/4" H	228.6 x 10.16 x .64 cm	Adjustable range: 59 1/2 x 79 1/2"	151.13 x 201.93 cm
	AKTON polymer surface measures 190.5 x 8.2 x 0.6cm			

## Vac-Pac® Covers




This pad is designed to line the interior of the Olympic Vac-Pac® (bean bag) positioner to help mitigate patient pressure and shear while in the lateral position. It is applicable for orthopedic, thoracic or urology procedures that dictate the patient be stable in the lateral position.

*Vac-Pac® is a registered trademark of Natus Medical Incorporated*

<b>40612</b>	Vac-Pac® Cover Medium Size 30	37 L x 28 W x 1/2" H	93.98 x 71.12 x 1.27 cm
<b>40612-31</b>	Vac-Pac® Cover Large Size 31	35 1/2 L x 38 3/4 W x 1/2" H	90.17 x 98.43 x 1.27 cm
<b>40613</b>	Vac-Pac® Cover T-Shape Size 35	15 L x 34 W x 1/2" H (bottom) 9 L x 23 W x 1/2" H (top)	38.1 x 86.36 x 1.27 cm (bottom) 22.86 x 58.42 x 1.27 cm (top)

## Breast Biopsy Table Pad

Improve patient comfort with this soft, polymer pad that helps protect skin and bony prominences while in the prone position during breast biopsy procedures.



**40302** 30 L x 22 3/4 W x 3/8" H with 9 3/8" D cutout 76.2 x 57.79 x .95 cm with 23.81 cm D cutout

### Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073  
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products  
feature our proprietary  
Akton® polymer



MADE IN  
USA

