

ACTION[®] ARM / SHOULDER / ELBOW POSITIONING



Hand / Wrist Support Pad

Reusable, polymer pad offers stable and comfortable hand and wrist positioning and replaces rolled up towels. Perfect for Carpal Tunnel, Arterial Cannulation, Dialysis and IV or A-line insertion support.

40321 11 L x 3 1/2 W x 2 1/4 H" 27.94 x 8.89 x 5.7 cm



Shown with Stabilizer Strap



Shown with Stabilizer Strap



40307 Small 40310 Medium 40306 Large

Contoured Armboard Pads

The contoured armboard design cradles patient's arm while providing full-arm protection of the bicep, elbow, and forearm. Gently raised sides to help keep arm securely positioned.

				3/4" Deep Middle Channel
40307	Small	16 L x 5 3/4 W x 2 1/4 H"	40.64 x 14.61 x 5.72 cm	
40310	Medium	20 L x 5 3/4 W x 2 1/4 H"	50.8 x 14.6 x 5.72 cm	
40306	Large	24 L x 5 3/4 W x 2 1/4 H"	60.96 x 14.6 x 5.72 cm	



Composite Armboard Pads

Three styles of armboard pads are available (Flat, Contoured, and Incline) complementing our full line of O.R. Replacement Table Pads. The five varying construction types offer a range of heights, pressure management options, and price points to meet every facility's needs.



► See our Support Surfaces Brochure for more details.



Arm / Hand Stabilizer Strap

This reusable strap features a soft polymer section to protect patient's skin while stabilizing arm/ hand. Strap is secured with attached hook and loop.

40315	Short	25 L x 2 W x 1/4 H" 12 1/2 L x 2 W x 1/4 H" (polymer area only)	63.5 x 5.08 x .64 cm
40317	Long	36 L x 2 W x 1/4 H" 15 1/2 L x 2 W x 1/4 H" (polymer area only)	91.44 x 5.08 x .64 cm



Shown with Stabilizer Strap



Armboard Pads

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas.

40300	Short	15 L x 4 W x 1/4 H"	38.1 x 10.16 x .64 cm
40301	Long	24 L x 4 W x 3/8 H"	60.96 x 10.16 x .95 cm
40309	Wide Short	15 L x 6 1/2 W x 1/4 H"	38.1 x 16.51 x .64 cm
40308	Wide Long	24 L x 6 1/2 W x 3/8 H"	60.96 x 16.51 x .95 cm



Ulnar / Forearm Protector

This pad helps to protect the patient's ulnar nerve and full forearm from damage due to pressure or shear. Designed to be wrapped around the arm and secured with the attached hook & loop.



Hook and loop attachment

40303 18 L x 13 W x 1/4 H" 45.72 x 33.02 x .64 cm

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products
feature our proprietary
Akton[®] polymer



Toboggan Pad

This pad is for use with the arm/ leg guard shield (toboggan). Together they secure the arm or leg position on the table and protect the patient from being leaned on or pinched during the procedure.

The polymer pad cushions and helps protect the arm's ulnar and radial nerves.

The small film pocket slips on to the shield and pad is secured by hook & loop attachment (included).

The shield is only sold with the set.



Shown with Heel Support



Hook and loop attachment



40130	Toboggan Pad	16 1/4 L x 8 1/4 W x 3/8 H"	41.27 x 20.96 x .95 cm
S40130*	Toboggan Pad Set	16 1/4 L x 8 1/4 W x 3/8 H"	41.27 x 20.96 x .95 cm
	Shield size =	15 1/2 L x 18 1/4 W"	39.37 x 46.4 cm

*Warranty covers pad only.

*Toboggan Armguard Shield is not available for sale or shipment outside of the United States and Canada.

Mayo Stand Cover

Use to prevent pressure and nerve damage when limb is positioned on the Mayo Stand during procedures.

5/8" polymer pad is sealed along three sides to a thick, durable film creating a pocket that slides easily on to the top of the Mayo stand.




40717 21 L x 18 1/2 W x 5/8 H" 53.3 x 47 x 1.6 cm

Hourglass Hand Table Pad

Reusable 1/2" polymer overlay designed for surgical procedures involving the forearm, elbow, or hand utilizing the hourglass hand table. Good for lengthy hand microsurgies requiring additional pressure and shear relieving support. Fits most standard hourglass hand tables.



40330 32 L x 24 W x 1/2 H" 81.28 x 60.96 x 1.27 cm
16"W middle section 40.6 cm middle section



AKTON polymer with soft, flexible film surface



High-density Foam

Hand Table Composite Pad

Designed for surgical procedures involving the forearm, elbow, or hand. Good for lengthy hand microsurgies requiring additional pressure and shear relieving support. Pad combines an AKTON polymer patient surface atop a high-density foam base featuring fully-welded, non-stitch construction. This pad is designed/constructed with the same pressure management benefits as the Action Products' "Original" O.R. table pads.

40320 34 L x 16 W x 2 3/8 H" 86.36 x 40.64 x 6.03 cm



Ulnar Nerve Protector

This pad helps to protect the patient's ulnar nerve while providing greater visibility and access for IV or A-line for anesthesiologist. Pad is placed around the elbow and secured with attached hook & loop. Maximum pad circumference is approximately 15.75" on the small and 19" on the large.

40304	Small	15 3/4 L x 6 W x 1/2 H"	40 x 15.24 x 1.27 cm
40304L	Large	19 L x 8 1/2 W x 1/2 H"	48.26 x 21.6 x 1.27 cm
40304P	Pediatric	10 L x 3 W x 1/4 H"	25.4 x 7.6 x .6 cm

