


ACTION® DISPOSABLE COVERS

- ✓ Soft, flexible film cover
- ✓ Fluid resistant
- ✓ Helps reduce clean up time
- ✓ Extends the life of the overlay



Disposable Overlay Covers



Action's fluid-resistant, non-sterile Disposable Overlay Covers are made of a soft, flexible film and are available in two styles. They are designed to help reduce exposure to fluids and to help extend the life of the overlay.

The "pillowcase" style is closed at one end and fully encases the overlay like a sheath.

The "fitted sheet" style covers the top of the pad and is held in place by two pockets at each end of the cover. This cover has a 4" film drape that runs the length of the cover on both sides.

Purchase in cases of 25.

Fitted Sheet Style		
40800	Made for overlay model 40100- medium overlay	Case of 25
40801	Made for overlay model 40101- full length overlay	Case of 25
40803	Made for overlay model 40103- torso with cut out	Case of 25
40805	Made for small overlay	Case of 25
Pillowcase Style		
40102	Made for table pad models 40100 & 40101	Case of 25



Disposable Chest Roll Covers

Action's fluid-resistant, non-sterile Disposable Chest Roll Covers are made of a soft, flexible film designed to extend the life of the chest roll and reduce exposure to fluids during use.

The slip-on covers mimic size and shape of the chest roll and have an overlapping film flap along the flat bottom.

Purchase in cases of 25.



Fits Action® Chest Roll		
40616DC	40616	Case of 25
40617DC	40617	Case of 25
40618DC	40618	Case of 25
40619DC	40619	Case of 25
40620DC	40620	Case of 25
40621DC	40621	Case of 25

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.



MADE IN USA



M3014-002