

# ACTION<sup>®</sup> PEDIATRIC POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions for even the smallest patients.

**Pressure injury protection, from head to heels.**

- ▶ Over 80% of the top U.S. pediatric hospitals use Action<sup>®</sup> products to protect their patients.

Neonates to Teens



## Head Positioners



### Donut Head Pads

Designed to protect and cradle the patient's head, these pads afford easy access for surgical area of the head and neck. It can be used for a variety of procedures including ENT, plastic, general and ophthalmic.

40202	Pediatric	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40219	Pediatric/ Neo-Natal	4 OD x 2 ID x 1"	10.2 OD x 6.4 ID x 2.5 cm
40203	Neo-Natal	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm
40202D	Pediatric with Center Dish	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 6.4 ID x 3.2 cm
40203D	Neo-Natal with Center Dish	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm

Teen / Adult sizes available



### Horseshoe Head Pads

The horseshoe configuration provides a clear air passageway and aids the anesthesiologist while patients are in a lateral and prone position, head to side.

40205	Pediatric	5 1/2 OD x 2 1/4 ID x 1 1/4"	14 OD x 5.7 ID x 3.2 cm
40220	Pediatric/ Neo-Natal	4 1/2 OD x 2 ID x 1"	11.4 OD x 5.1 ID x 2.5 cm
40206	Neo-Natal	3 1/4 OD x 1 1/2 ID x 3/4"	8.3 OD x 3.8 ID x 2 cm

Teen / Adult sizes available

### Head Pad with Center Dish

Designed for smaller patients, this moderately dished, square head pad may be used to minimize pressure on the occiput, neck, head or ears of the patient as needed.

40221 Pediatric 6 x 6 x 1" 15.2 x 15.2 x 2.5 cm

Teen / Adult sizes available



### Ophthalmic Cradle Headrest

Designed originally for cataract surgery, this versatile pad helps prevent risk of alopecia - especially during long surgeries. Cradling of patient's head provides stability in addition to easy surgeon/nurse access. It is ideal for many procedures including ENT, plastic surgery and many ophthalmic cases.

40217P Pediatric 7 1/2 x 8 1/2 x 2 1/4" 19 x 21.6 x 5.7 cm

Teen / Adult sizes available



### Forehead Protector

The Forehead Protector is semi-circular with a flat bottom. It is used to protect the infant or child patient's forehead from pressure exerted by oxygen tubes.

40208 2 1/2 x 1 5/8 x 5/8" 6.3 x 4 x 1.6 cm



From neonates to teens, Action Products has a full range of pediatric positioners. If we don't have the size you need we can fabricate pads to fit your specifications.

Email [service@actionproducts.com](mailto:service@actionproducts.com)  
Call 800.228.7763 or 301.797.1414

**Action Products, Inc.**

800.228.7763 • 301.797.1414 • Fax 877.733.2073  
[www.actionproducts.com](http://www.actionproducts.com) • [service@actionproducts.com](mailto:service@actionproducts.com)

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.

ACTION products  
feature our proprietary  
Akton<sup>®</sup> polymer



## Body Positioners



### Ulnar Nerve Protector

Provides pressure injury protection for pediatric elbows and ulnar nerves. Pad has a hook and loop attachment for adjustable, secure fit.

<b>40304P</b>	Pediatric	10 x 3 x ¼"	25.4 x 7.6 x .6 cm
Teen / Adult sizes available			



### Universal Spot Protectors

This series of thin, round pads is designed to provide adaptive use padding options for those hard-to-protect areas unique to each patient.

Thin profile, 1/4" thick polymer



<b>40405</b>	Mini	1 ½" D x ¼"	4.1 D x .6 cm	Pkg of 12
<b>40406</b>	Small	3 ¼" D x ¼"	8.3 D x .6 cm	Pkg of 12
<b>40407</b>	Large	5 ½" D x ¼"	14 D x .6 cm	Pkg of 6



Bassinet Pad



### NICU Pads

Provides pressure and shear management surface for pediatric and neo-natal patients.

<b>40120L</b>	Isolette Pad	26 x 16 x ½"	66 x 40.6 x 1.3 cm
<b>40121</b>	Isolette Pad	29 x 19 x ½"	74 x 48.2 x 1.3 cm
<b>40120</b>	Bassinet Pad	26 x 12 x ½"	66 x 30.5 x 1.3 cm

## Pediatric Positioner Set

This versatile 4-piece set is designed to position and protect newborns, infants and toddlers during a variety of procedures. Placed flat and used alone, the channeled positioner can help stabilize patients while in supine. Used in combination with a chest roll under the pad at the thoracic area, the chest is extended for anterior chest procedures. Additionally, when used with the small chest roll under the channeled positioner at the thorax and large chest roll at the sacrum, access for posterior supine cases such as spinal tumor resection is enhanced.

Visit our website for more information on this set.



**40614** Four piece pediatric positioning set  
Includes: 40205, 40642, 40646, 40647

**40647** Rectangular Channel Pad  
20 x 9 x 2" 50.8 x 22.9 x 5 cm



## Pediatric Chest Rolls

Designed for smaller patients, this line of flat bottom chest rolls provide stable support in the prone position. They are available in a wide range of heights and lengths to accommodate virtually any patient.

<b>40639</b>	6 x 2 x 2"	15 x 5 x 5 cm
<b>40640</b>	8 x 2 x 2"	20 x 5 x 5 cm
<b>40641</b>	10 x 2 x 2"	25.5 x 5 x 5 cm
<b>40642</b>	12 x 2 x 2"	30.5 x 5 x 5 cm
<b>40643</b>	6 x 3 x 3"	15 x 7.5 x 7.5 cm
<b>40644</b>	8 x 3 x 3"	20 x 7.5 x 7.5 cm
<b>40645</b>	10 x 3 x 3"	25.5 x 7.5 x 7.5 cm
<b>40646</b>	12 x 3 x 3"	30.5 x 7.5 x 7.5 cm

Teen / Adult sizes available



## Pediatric/Neo-Natal Positioners

These positioners are used to position infants without rolling wash-cloths and other items.

<b>40648</b>	6 x 1 ½ x ¾"	15 x 3.8 x 2 cm
<b>40649</b>	6 x 1 x 1"	15 x 2.5 x 2.5 cm
<b>40650</b>	6 x 1 ½ x 1 ½"	15 x 3.8 x 3.8 cm

Teen / Adult sizes available



## Patient Axillary Positioner

AORN encourages the use of Axillary Rolls in place of sandbags or IV bags. These extremely versatile positioners are good for many procedures.

<b>40602</b>	Small	10 x 4 ¼ x 2"	25.4 x 11.4 x 4 cm
--------------	-------	---------------	--------------------

Teen / Adult sizes available