

ACTION[®] SUPINE POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions such as supine.
Pressure injury protection, from head to heels.

Head Positioners

► Please see our **Pediatric Product Sheet** for our extensive range of Pediatric Positioning Products.



40201

40201D

40218

Donut Head Pads

Designed to protect and cradle the patient's head, these pads afford easy access for surgical area of the head and neck. It can be used for a variety of procedures including ENT, plastic, general and ophthalmic.

40201	Adult	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40201D	Adult w/ Center Dish	8 OD x 3 ID x 1 3/4"	20.5 OD x 7.5 ID x 4.5 cm
40218	Adult Wide	8 OD x 4 1/2 ID x 1 3/4"	20.5 OD x 11.4 ID x 4.5 cm

Supine Head Pad

The Action[®] Supine Head Pad may be used for most supine positions. After placing patient's head in the pad, the centering dish conforms to the patient's head for support and pressure prevention.



40200	9 x 10 x 2 1/4"	23 x 25 x 5.7 cm
--------------	-----------------	------------------

Head Pad with Center Dish

This moderately dished, square head pad may be used to minimize pressure on the occiput, neck or head of the patient as needed.



40207	9 x 10 x 1"	23 x 25 x 2.5 cm
--------------	-------------	------------------



Contoured Head Pad

The pads contoured design provides optimum stability of the patient's head while in the supine position.

40213	7 x 7 x 3 1/4" (1 1/2" contoured)	18 x 18 x 8 cm
--------------	-----------------------------------	----------------

Ophthalmic Headrests

The Action[®] Ophthalmic Headrest Pad includes a contoured cutout for stable head placement. The robust edges are designed to support the surgeon's hands if needed during surgery.



40209	Large	12 x 10 x 3 7/8"	31 x 25.4 x 10 cm
40210	Medium	12 x 10 x 2 7/8"	31 x 25.4 x 7.3 cm

Ophthalmic Cradle Headrest

Designed originally for cataract surgery, this versatile pad helps prevent risk of alopecia, especially during long surgeries.



40217	Adult	9 x 10 1/2 x 2 5/8"	23 x 26.7 x 6.7 cm
--------------	-------	---------------------	--------------------

Forehead Protector

The Forehead Protector is semi-circular with a flat bottom. It is used to protect the infant, child or adult patient's forehead from pressure exerted by oxygen tubes.



40208	2 1/2 x 1 5/8 x 5/8"	6.3 x 4 x 1.6 cm
--------------	----------------------	------------------

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.

ACTION products
feature our proprietary
Akton[®] polymer



Arm/Elbow/Hand Positioners

Ulnar Forearm Protector

Designed to be wrapped around the arm and secured with the Hook & Loop attachment as the patient's arm is tucked to their side.



40303 18 x 13 x ¼" 46 x 33 x .6 cm

Ulnar Nerve Protector

Designed to allow a greater view of the arm while permitting easy access by the anesthesiologist. Pad is placed over elbow and securely attached with hook & loop attachment.



40304 Small 15 7/8 x 6 x ½" 40.3 x 15 x 1.3 cm

40304L Large 19 x 8 ½ x ½" 48 x 21.6 x 1.3 cm

Hand/Wrist Support

Superior hand and wrist stabilization and protection. Perfect for Carpal Tunnel, Arterial Cannulation and Dialysis and other procedures.




40321 Hand/Wrist Support 11 x 3 ½ x 2 ¼" 28 x 9 x 5.7 cm

40315 Strap 25 x 2 x ¼" 63.5 x 5 x .635 cm

40317 Strap 36 x 2 x ¼" 91.4 x 5 x .635 cm

Toboggan Pad

This pad is intended for use in conjunction with arm guard shield (toboggan). Designed to prevent patient's arm or leg from extending over the side of the table and to prevent leaning on or against the patient's extremity. *The shield is only sold with the set.*




40130 Toboggan Pad 16 ¼ x 8 ¼ x ¾" 41.3 x 21 x 1 cm

S40130* Toboggan Pad Set 16 ¼ x 8 ¼ x ¾" 41.3 x 21 x 1 cm

Shield size = 15 ¼ x 18 ¼" 38.7 x 46.4 cm

*Warranty covers pad only.




Contoured Armboard Pads

These pads offer excellent full-arm coverage. Contoured design cradles patient's arm, protects bicep, elbows and forearms.

40307 Small 16 x 5 ¾ x 2 ¼" 41 x 14.6 x 5.7 cm

40310 Medium 20 x 5 ¾ x 2 ¼" 51 x 14.6 x 5.7 cm

40306 Large 24 x 5 ¾ x 2 ¼" 61 x 14.6 x 5.7 cm



Armboard Pads

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas.

40300 Short 15 x 4 x ¼" 38 x 10.2 x .6 cm

40301 Long 24 x 4 x ¾" 61 x 10.2 x 1 cm

40309 Wide Short 15 x 6 ½ x ¼" 38 x 16.5 x .6 cm

40308 Wide Long 24 x 6 ½ x ¾" 61 x 16.5 x 1 cm

40315 Strap 25 x 2 x ¼" 63.5 x 5 x .635 cm

40317 Strap 36 x 2 x ¼" 91.4 x 5 x .635 cm



Mayo Stand Cover

Use to prevent pressure and nerve damage when limb is placed on the Mayo Stand for positioning during orthopedic procedures.



40717 21 x 18 ½ x ¾" 53.3 x 47 x 1.6 cm

Leg/Ankle/Heel Positioners

Dome-Shaped Positioners

40603	14 x 6 x 3"	36 x 15 x 7.6 cm
40603L	20 x 6 x 3"	50.8 x 15.2 x 7.6 cm



Semi-circular design with flat bottom configuration provides pressure protection under the knees while the patient is in a supine position. Helps to protect patient against foot drop for patients in prone position. Can be used under the metatarsal or under the foot, against the board, or in the semi-fowlers position.

Trapezoid Positioner

Designed for surgical procedures in which the patient must have clear leg separation.

40630	10 x 7 x 3½"	25.4 x 17.8 x 8.9 cm
--------------	--------------	----------------------



Positioning Strap

Welded seams and 4-way stretch microbe-resistant fabric make this a must have for any operating room. Radiolucent and MRI safe. Hook and loop attachments secure the strap to standard tables.

This strap is not meant to restrain patient.

40713	57 x 4 x ¼"	144.8 x 10.2 x .6 cm	Adjustable range: 28 x 43"	71.12 x 109.22 cm
40713L	90 x 4 x ¼"	228.6 x 10.2 x .6 cm	Adjustable range: 59½ x 79½"	151.13 x 201.93 cm



Heel Support

40502	6 x 3 ¼" x 1 ¾"	15 x 8.3 x 4.5 cm
--------------	-----------------	-------------------



Place patient's heel in cup cut-out. Slit in back of cup secures the patient's Achilles tendon area.



Heel Protector Pad

40401	12 ⅞ x 7 ⅞ x ½"	32.7 x 20 x 1.3 cm
--------------	-----------------	--------------------



Pads should be used for any procedure that is lengthy and involves pressure on the patient's heel area.



Heel Comfort Cup Set

40520 (Set of 2)	4 ¼ x 3 5/8 x ½"	10.8 x 9.2 x 1.3 cm
-------------------------	------------------	---------------------



Provides comfort and support for heels. May be used to complement anti-embolism stockings. Appropriate for recovery and radiology. Comes complete with polymer fastening strap. Universal sizing.

Advantages of AKTON® Polymer

- Clinically proven
- Distributes weight evenly
- Distributes heat evenly
- Will not leak or bottom out
- Will not support bacterial growth
- Imaging compatible
- Fast and easy cleaning
- Reusable

- Cost Effective
- Tissue Equivalent
- Hypo-Allergenic
- Free of plasticizers, silicone, DEHP
- Not made with natural rubber latex





SUPINE POSITIONING PRODUCTS

Support Surfaces & Overlays

Support Surfaces

Benefits you can expect from every Action® support surface:

- Imaging compatible
- Meets ergonomic weight limits for hospital staff
- Not made with natural rubber latex, silicone and plasticizer free
- Static resistive
- Easy to clean
- 100% no-stitch construction
- Fluid/ waterproof with fully sealed seams
- Available in multiple thicknesses
- Coordinating armboards available
- Table Pads for all major Brands



Overlays

The Action® O.R. Overlay is a versatile pad. Please visit www.actionproducts.com to see our full line of Action® Overlays.

40105	Small	20 x 20 x 1/2"	50.8 x 50.8 x 1.3 cm
40100	Medium	46 x 20 x 1/2"	117 x 50 x 1.3 cm
40109	Medium Large	63 x 20 x 1/2"	160 x 50 x 1.3 cm
40101	Full Length	72 x 20 x 1/2"	183 x 50 x 1.3 cm



Call customer service
800.228.7763 or 301.797.1414,
or email service@actionproducts.com

Performance Series

Our most popular table pad for high turnover Operating Rooms

- 4-Way Stretch Fabric Cover
- Use for heavy frames
- Encased 3/8" head-to-toe layer of AKTON viscoelastic polymer provides immersion benefits by distributing weight
- Bacteria / Microbe resistant solution in fibers of fabric

Performance Plus Series

- Exposed continuous layer of our proprietary AKTON polymer surface provides excellent shear benefits and will not support bacterial growth
- Soft viscofoam mid-layer enhances patient envelopment and comfort
- Great choice for extended cases

Angel Series

- Our most budget friendly alternative to the "basic" pad
- Upgraded viscofoam/ high density foam construction
- Use with Action® Overlays
- Comfort for non-anesthetized patient procedures

► See our **Support Surface Brochure** for more detailed information on our full range of table pads.

Patient Protection for Supine Procedures including:

Abdominoplasty	Cholecystectomy	Septoplasty
Blepharoplasty	Hernia Repair	Tonsillectomy
Breast Biopsy	Knee Arthroscopy	
Cardiovascular	Rhytidectomies	



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.



MADE IN
USA



M3022-001