

# ACTION<sup>®</sup> PRONE POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions such as prone. **Pressure injury protection, from head to heels.**



## Head Positioners & Overlays



### Prone Headrests

Our popular Prone Headrest is designed to help keep the patient's head in a straight, neutral position during prone procedures.

<b>40216M</b>	Medium	8 ¾ x 7 ¾ x 4 ¾"	22.2 x 19.7 x 12 cm
<b>40216L</b>	Large	10 ½ x 9 x 5 ¼"	27 x 22.0 x 13.4 cm

### Horseshoe Head Pad

The horseshoe configuration provides a clear air passageway and aids the anesthesiologist while patients are in a lateral and prone position, head to side.

<b>40204</b>	Adult	8 OD x 3 ID x 1 ¾"	20.5 OD x 7.5 ID x 4.5 cm
<b>40211</b>	Adult, High	8 OD x 3 ID x 3"	20.5 OD x 7.5 ID x 7.5 cm

### Overlays

The Action<sup>®</sup> O.R. Overlay is a versatile pad used for a wide range of procedures.

<b>40105</b>	Small	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm
<b>40100</b>	Medium	46 x 20 x ½"	117 x 50 x 1.3 cm
<b>40101</b>	Full Length	72 x 20 x ½"	183 x 50 x 1.3 cm
<b>40103</b>	Perineal Cutout	46 x 20 x ½" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
<b>40104</b>	3-piece Set	See A, B and C	
<b>40104A</b>	Head	9 x 20 x ½"	22.9 x 50.8 x 1.3 cm
<b>40104B</b>	Torso, w/ perineal cutout	39 ½ x 20 x ½"	100.3 x 50.8 x 1.3 cm
<b>40104C</b>	Foot, w/ perineal cutout	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm

## Arm Positioners



### Contoured Armboard Pads

These pads offer excellent full-arm coverage. Contoured design cradles patient's arm, protects bicep, elbows and forearms.

<b>40307</b>	Small	16 x 5 ¾ x 2 ¼"	41 x 14.6 x 5.7 cm
<b>40310</b>	Medium	20 x 5 ¾ x 2 ¼"	51 x 14.6 x 5.7 cm
<b>40306</b>	Large	24 x 5 ¾ x 2 ¼"	61 x 14.6 x 5.7 cm



### Armboard Pads

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas.

<b>40317</b>	Strap	36 x 2 x ½"	91.4 x 5 x 1 cm
<b>40315</b>	Strap	25 x 2 x ¼"	63.5 x 5 x .635 cm
<b>40300</b>	Short	15 x 4 x ¼"	38 x 10.2 x .6 cm
<b>40301</b>	Long	24 x 4 x ¾"	61 x 10.2 x 1 cm
<b>40309</b>	Wide Short	15 x 6 ½ x ¼"	38 x 16.5 x .6 cm
<b>40308</b>	Wide Long	24 x 6 ½ x ¾"	61 x 16.5 x 1 cm

## Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073  
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2021 Action Products, Inc.

ACTION products feature our exclusive Akton<sup>®</sup> polymer



## Patient protection for the following procedures:

- |                        |  |
|------------------------|--|
| Achilles Tendon Repair | Posterior Cervical Fusion              |
| Hemorrhoidectomy       | Posterior Varicose Vein Avulsion       |
| Lipoma Excision        | Repair of Haglund's Deformity          |
| Lumbar Laminectomy     | Suction Assisted Lipectomy of Buttocks |



## Body Positioners



### Contoured Chest Rolls

These flat bottom chest rolls are contoured to conform to an adult's rib cage. Pad provides maximum stability and great utility. Available in three different sizes.

<b>40622</b>	Small	12 x 6 ¼ x 5"	30.5 x 16 x 13 cm
<b>40623</b>	Medium	16 x 6 ¼ x 5"	40.6 x 16 x 13 cm
<b>40624</b>	Large	20 x 6 ¼ x 5"	50.8 x 16 x 13 cm



### Chest Rolls

This line of flat bottom chest rolls is produced in a wide range of sizes to accommodate virtually any person.

Adult			Pediatric	
<b>40616</b>	12 ½ x 4 ¾ x 3 ¾"	31.8 x 11.1 x 8.6 cm	<b>40639</b>	6 x 2 x 2" / 15 x 5 x 5 cm
<b>40617</b>	16 ½ x 4 ¾ x 3 ¾"	41.9 x 11.1 x 8.6 cm	<b>40640</b>	8 x 2 x 2" / 20 x 5 x 5 cm
<b>40618</b>	20 ½ x 4 ¾ x 3 ¾"	52.1 x 11.1 x 8.6 cm	<b>40641</b>	10 x 2 x 2" / 25.5 x 5 x 5 cm
<b>40619</b>	12 x 6 x 5"	30.5 x 15 x 13 cm	<b>40642</b>	12 x 2 x 2" / 30.5 x 5 x 5 cm
<b>40620</b>	16 x 6 x 5"	40.6 x 15 x 13 cm	<b>40643</b>	6 x 3 x 3" / 15 x 7.5 x 7.5 cm
<b>40621</b>	20 x 6 x 5"	50.8 x 15 x 13 cm	<b>40644</b>	8 x 3 x 3" / 20 x 7.5 x 7.5 cm
			<b>40645</b>	10 x 3 x 3" / 25.5 x 7.5 x 7.5 cm
			<b>40646</b>	12 x 3 x 3" / 30.5 x 7.5 x 7.5 cm



### Chest Positioners

These positioners have a contoured design patterned to uniformly conform to a patient's shoulder-to-chest area.

<b>40625</b>	Left	28 x 23 ¾ x 8 x 3 x 2 ½"	71.1 x 59.4 x 20.3 x 7.6 x 6.4 cm
<b>40626</b>	Right	28 x 23 ¾ x 8 x 3 x 2 ½"	71.1 x 59.4 x 20.3 x 7.6 x 6.4 cm

### Dome-Shaped Positioners

Semi-circular design with flat bottom configuration provides pressure protection under the knees while the patient is in a supine position. Helps to protect patient against foot drop for patients in prone position. Can be used under the metatarsal or under the foot, against the board, or in the semi-fowlers position.



<b>40603</b>	14 x 6 x 3"	36 x 15 x 7.6 cm
<b>40603L</b>	20 x 6 x 3"	50.8 x 15.2 x 7.6 cm