

ACTION[®] PRONE POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions such as prone.
Pressure injury protection, from head to heels.

Head Positioners & Overlays



Prone Headrests

Our popular Prone Headrest is designed to support the patient's head in a straight, neutral position during prone procedures. Sculpted central cutout is designed to avoid pressure around eyes, nose, and face. Soft, Akton[®] polymer surface helps reduce shear at the forehead and chin support area. Channels present in top and base of pad provide tube access for the anesthesiologist.


40216M	Medium	8 ¾ x 7 ¾ x 4 ¾"	22.2 x 19.7 x 12 cm
40216L	Large	10 ¾ x 9 x 5 ¼"	27 x 22.0 x 13.4 cm

Overlays

The Action[®] O.R. Overlay is a versatile pad used for a wide range of procedures.

40105	Small	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm
40100	Medium	46 x 20 x ½"	117 x 50 x 1.3 cm
40109	Medium Large	63 x 20 x ½"	160 x 50 x 1.3 cm
40101	Full Length	72 x 20 x ½"	183 x 50 x 1.3 cm
40103	Perineal Cutout	46 x 20 x ½" w/ cutout	117 x 50.8 x 1.3 cm w/ cutout
40104	3-piece Set	See A, B and C	
40104A	Head	9 x 20 x ½"	22.9 x 50.8 x 1.3 cm
40104B	Torso, w/ perineal cutout	39 ½ x 20 x ½"	100.3 x 50.8 x 1.3 cm
40104C	Foot, w/ perineal cutout	20 x 20 x ½"	50.8 x 50.8 x 1.3 cm

Arm Positioners



Contoured Armboard Pads

These pads offer excellent full-arm coverage. Contoured design cradles patient's arm, protects bicep, elbows and forearms.


40307	Small	16 x 5 ¾ x 2 ¼"	41 x 14.6 x 5.7 cm
40310	Medium	20 x 5 ¾ x 2 ¼"	51 x 14.6 x 5.7 cm
40306	Large	24 x 5 ¾ x 2 ¼"	61 x 14.6 x 5.7 cm



Armboard Pads

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas.

40300	Short	15 x 4 x ¼"	38 x 10.2 x .6 cm
40301	Long	24 x 4 x ¾"	61 x 10.2 x 1 cm
40309	Wide Short	15 x 6 ½ x ¼"	38 x 16.5 x .6 cm
40308	Wide Long	24 x 6 ½ x ¾"	61 x 16.5 x 1 cm
40315	Strap	25 x 2 x ¼"	63.5 x 5 x .635 cm
40317	Strap	36 x 2 x ¼"	91.4 x 5 x .635 cm



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products
feature our proprietary
Akton[®] polymer





PRONE POSITIONING PRODUCTS

Patient protection for the following procedures:

Achilles Tendon Repair	Posterior Cervical Fusion
Hemorrhoidectomy	Posterior Varicose Vein Avulsion
Lipoma Excision	Repair of Haglund's Deformity
Lumbar Laminectomy	Suction Assisted Lipectomy of Buttocks



Body Positioners



40616 40617 40618 40621 40620 40619

Chest Rolls


This line of flat bottom chest rolls is designed to provide stable support in the prone position. They are available in a wide range of heights and lengths to accommodate virtually any patient, adult to neonate.



40643 40644 40645 40646 40639 40640 40641 40642

Adult		
40616	12 1/2 L x 4 3/4 W x 3 3/8 H	31.75 x 11.11 x 8.57 cm
40617	16 1/2 L x 4 3/4 W x 3 3/8 H	41.91 x 11.11 x 8.57 cm
40618	20 1/2 L x 4 3/4 W x 3 3/8 H	52.07 x 11.11 x 8.57 cm
40619	12 L x 6 W x 5 H	30.48 x 15.24 x 12.7 cm
40620	16 L x 6 W x 5 H	40.64 x 15.24 x 12.7 cm
40621	20 L x 6 W x 5 H	50.8 x 15.24 x 12.7 cm

Pediatric		
40639	6 L x 2 W x 2 H	15.24 x 5.08 x 5.08 cm
40640	8 L x 2 W x 2 H	20.32 x 5.08 x 5.08 cm
40641	10 L x 2 W x 2 H	25.4 x 5.08 x 5.08 cm
40642	12 L x 2 W x 2 H	30.48 x 5.08 x 5.08 cm
40643	6 L x 3 W x 3 H	15.24 x 7.62 x 7.62 cm
40644	8 L x 3 W x 3 H	20.32 x 7.62 x 7.62 cm
40645	10 L x 3 W x 3 H	25.4 x 7.62 x 7.62 cm
40646	12 L x 3 W x 3 H	30.48 x 7.62 x 7.62 cm



Chest Positioners

These positioners have a contoured design patterned to uniformly conform to a patient's shoulder to chest area, with shaped edges for the axilla. They provide upper body positioning while alleviating pressure on the thorax area.

40625	Left	28 L x 23 3/8 L x 8 x 3 x 2 1/2 H	71.1 L x 59.4 L x 20.3 x 7.6 x 6.4 H cm
40626	Right	28 L x 23 3/8 L x 8 x 3 x 2 1/2 H	71.1 L x 59.4 L x 20.3 x 7.6 x 6.4 H cm

Dome-Shaped Positioners

Semi-circular design with flat bottom configuration provides pressure protection under the knees while the patient is in a supine position. Helps to protect patient against foot drop for patients in prone position. Can be used under the metatarsal or under the foot, against the board, or in the semi-fowlers position.



40603	14 x 6 x 3	36 x 15 x 7.6 cm
40603L	20 x 6 x 3	50.8 x 15.2 x 7.6 cm

Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073

www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.



MADE IN USA



M3024-002