



LATERAL POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions such as lateral.

Pressure injury protection, from head to heels.


Head Positioners




Head Pad with Center Dish

This moderately dished head pad reduces pressure on the side of the patient's face and ear.

40207	Adult	9 x 10 x 1"	23 x 25 x 2.5 cm
40221	Pediatric	6 x 6 x 1"	15.2 x 15.2 x 2.5 cm




Arm/Elbow/Hand Positioners



Contoured Armboard Pads

These pads offer excellent full-arm coverage. Contoured design cradles patient's arm, protects bicep, elbows and forearms.

40307	Small	16 x 5 3/4 x 2 1/4"	41 x 14.6 x 5.7 cm
40310	Medium	20 x 5 3/4 x 2 1/4"	51 x 14.6 x 5.7 cm
40306	Large	24 x 5 3/4 x 2 1/4"	61 x 14.6 x 5.7 cm



Armboard Pads & Stabilizer Straps

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas. The polymer Stabilizer Strap secures and protects the arm.

40300	Short	15 x 4 x 1/4"	38 x 10.2 x .6 cm
40301	Long	24 x 4 x 3/8"	61 x 10.2 x 1 cm
40309	Wide Short	15 x 6 1/2 x 1/4"	38 x 16.5 x .6 cm
40308	Wide Long	24 x 6 1/2 x 3/8"	61 x 16.5 x 1 cm
40315	Strap	25 x 2 x 1/4"	63.5 x 5 x .635 cm
40317	Strap	36 x 2 x 1/4"	91.4 x 5 x .635 cm



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.

ACTION products
feature our proprietary
Akton® polymer





LATERAL POSITIONING PRODUCTS

Body Positioners



Patient Axillary Positioner

AORN encourages the use of Axillary Rolls in place of sandbags or IV bags. These extremely versatile positioners are good for many procedures.

40602	Small	10 x 4 ¼ x 2"	25.4 x 11.4 x 5 cm
40605	Large	14 x 4 x 2"	35.5 x 10.2 x 5 cm
40615	Extra Long	18 x 4 x 2"	46 x 10 x 5 cm



Leg/Ankle/Heel Positioners



Heel Protector Pad

Pads designed for any procedure that requires pressure protection on both the patient's heel and ankle. Hook and loop closure accommodates a variety of leg and ankle sizes.

40401	12 7/8 x 7 7/8 x 1/2"	32.7 x 20 x 1.3 cm
--------------	-----------------------	--------------------



Lateral Positioner

Accommodates both shoulder and upper arm. Provides pressure relief and nerve protection. Very suitable for chest, orthopedic or kidney procedures. Lightweight foam base with AKTON polymer pad overlay.

40610	28 x 19 1/2 x 6 - 1"	71.1 x 49.5 x 15.2 - 2.5 cm
--------------	----------------------	-----------------------------



Positioning Strap

Welded seams and waterproof film combined with a soft patient-facing surface make this a must have for the operating room. Radiolucent and MRI safe. Hook and loop attachments secure the strap to standard tables. *This strap is not meant to restrain patient.*

40713	57 x 4 x 1/4"	144.8 x 10.2 x .6 cm	Adjustable range: 28 x 43"	71.12 x 109.22 cm
40713L	90 x 4 x 1/4"	228.6 x 10.2 x .6 cm	Adjustable range: 59 1/2 x 79 1/2"	151.13 x 201.93 cm



Heel & Medical Elbow Protectors

Recommended for those clients "at risk" or diagnosed with compromised skin disorders. Full polymer strap offers enhanced, adjustable protection for heel or elbow.

20401	1/2" Thick
20401A	1/4" Thick



Advantages of AKTON® Polymer

- Clinically proven
- Distributes weight evenly
- Distributes heat evenly
- Will not leak or bottom out
- Will not support bacterial growth
- Imaging compatible
- Fast and easy cleaning
- Reusable

- Cost Effective
- Tissue Equivalent
- Hypo-Allergenic
- Free of plasticizers, silicone, DEHP
- Not made with natural rubber latex



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073

www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.





LATERAL POSITIONING PRODUCTS

Support Surfaces & Overlays

Support Surfaces

Benefits you can expect from every Action® support surface:

- Imaging compatible
- Meets ergonomic weight limits for hospital staff
- Not made with natural rubber latex, silicone and plasticizer free
- Static resistive
- Easy to clean
- 100% no-stitch construction
- Fluid/ waterproof with fully sealed seams
- Available in multiple thicknesses
- Coordinating armboards available
- Table Pads for all major Brands



Call customer service
800.228.7763 or 301.797.1414,
or email service@actionproducts.com



Overlays

The Action® O.R. Overlay is a versatile pad. Please visit www.actionproducts.com to see our full line of Action® Overlays.

40105	Small	20 x 20 x 1/2"	50.8 x 50.8 x 1.3 cm
40100	Medium	46 x 20 x 1/2"	117 x 50 x 1.3 cm
40109	Medium Large	63 x 20 x 1/2"	160 x 50 x 1.3 cm
40101	Full Length	72 x 20 x 1/2"	183 x 50 x 1.3 cm

Performance Series

Our most popular table pad for high turnover Operating Rooms

- 4-Way Stretch Fabric Cover
- Use for heavy frames
- Encased 3/8" head-to-toe layer of AKTON viscoelastic polymer provides immersion benefits by distributing weight
- Bacteria / Microbe resistant solution in fibers of fabric

Performance Plus Series

- Exposed continuous layer of our proprietary AKTON polymer surface provides excellent shear benefits and will not support bacterial growth
- Soft viscofoam mid-layer enhances patient envelopment and comfort
- Great choice for extended cases

Angel Series

- Our most budget friendly alternative to the "basic" pad
- Upgraded viscofoam/ high density foam construction
- Use with Action® Overlays
- Comfort for non-anesthetized patient procedures

► See our **Support Surface Brochure** for more detailed information on our full range of table pads.

Patient Protection for Lateral Procedures including:

- | | |
|------------------------|-----------------------------------|
| Colonoscopy | Hip Fracture |
| Intramedullary Rodding | Lumbar Epidural Steroid Injection |
| Lung Wedge Resection | Nephrectomy |
| Nephrolithotomy | Thoracoscopy |
| Tonsillectomy | |



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2023 Action Products, Inc.

ACTION products
feature our proprietary
Akton® polymer



M3026-001