



LATERAL POSITIONING PRODUCTS

Our reusable products are engineered to help preserve skin integrity and sustain clinically necessary positions such as lateral.

Pressure injury protection, from head to heels.


Head Positioners



Head Pad with Center Dish

This moderately dished head pad designed for lateral position reduces pressure on the side of the patient's face and ear.

40207	Adult	9 L x 10 W x 1 H"	22.9 x 25.4 x 2.5 cm
40221	Pediatric	6 L x 6 W x 1 H"	15.2 x 15.2 x 2.5 cm



Arm/Elbow/Hand Positioners




Contoured Armboard Pads

The contoured armboard design cradles patient's arm while providing full-arm protection of the bicep, elbow, and forearm. Gently raised sides to help keep arm securely positioned.

40307	Small	16 L x 5 3/4 W x 2 1/4 H"	40.64 x 14.61 x 5.72 cm
40310	Medium	20 L x 5 3/4 W x 2 1/4 H"	50.8 x 14.6 x 5.72 cm
40306	Large	24 L x 5 3/4 W x 2 1/4 H"	60.96 x 14.6 x 5.72 cm



3/4" Deep
Middle Channel



Armboard Pads & Stabilizer Straps

For use on standard armboard table extensions to protect the patient's bicep, elbow and forearm areas. The reusable Stabilizer Strap features a soft polymer section to protect patient's skin while stabilizing arm/ hand.

40300	Short	15 L x 4 W x 1/4 H"	38.1 x 10.16 x .64 cm
40301	Long	24 L x 4 W x 3/8 H"	60.96 x 10.16 x .95 cm
40309	Wide Short	15 L x 6 1/2 W x 1/4 H"	38.1 x 16.51 x .64 cm
40308	Wide Long	24 L x 6 1/2 W x 3/8 H"	60.96 x 16.51 x .95 cm
40315	Short Strap	25 L x 2 W x 1/4 H"	63.5 x 5.08 x .64 cm
40317	Long Strap	36 L x 2 W x 1/4 H"	91.44 x 5.08 x .64 cm



Action Products, Inc.

800.228.7763 • 301.797.1414 • Fax 877.733.2073
www.actionproducts.com • service@actionproducts.com

ACTION and AKTON are registered trademarks of Action Products, Inc. © 2024 Action Products, Inc.

ACTION products
feature our proprietary
Akton® polymer



Body Positioners



Lateral Positioner

Designed to protect the patient while in the lateral position. Accommodates both shoulder and upper arm. Provides pressure relief and nerve protection. Includes a lightweight foam base with AKTON polymer pad overlay.

40610 28 L x 19 ½ W x 6 - 1" H 71.12 x 49.53 x 15.24 - 2.54 cm

Patient Axillary Positioner

AORN encourages the use of Axillary Rolls in place of sandbags or IV bags. These extremely versatile positioners are good for many procedures.

40602	Small	10 L x 4 ¼ W x 1 ½" H	25.4 x 10.8 x 4.76 cm
40605	Long	14 L x 4 W x 2" H	35.56 x 10.16 x 5.08 cm
40615	Extra Long	18 L x 4 W x 2" H	45.72 x 10.16 x 5.08 cm



Positioning Strap

Welded seams and fluid-resistant film combined with a soft patient facing surface make this a must have for the operating room. Radiolucent and MRI safe. Hook and loop attachments secure the strap to standard tables. *This strap is not meant to restrain patient.*

40713	57 L x 4 W x ¼" H	144.78 x 10.16 x .64 cm	Adjustable range: 28 x 43"	71.12 x 109.22 cm
40713L	90 L x 4 W x ¼" H	228.6 x 10.16 x .64 cm	Adjustable range: 59½ x 79½"	151.13 x 201.93 cm

Leg/Ankle/Heel Positioners

Heel Protector Pad

This "clamshell" style pad provides full coverage ankle and heel protection during lengthy surgeries.

Interior surface of the pad is permanently labeled to ensure proper foot placement.

Easy to use and secure with adjustable hook and loop attachment. Accommodates approximately, a 15" ankle circumference.

Use alone in supine and lateral position, or combine with Ankle Strap Stirrup and Candy Cane pads for lithotomy position.



40401 12 ½ L x 7 ¾ W x ½ H" 32 x 20 x 1.2 cm



Vac-Pac® Covers

This pad is designed to line the interior of the Olympic Vac-Pac® (bean bag) positioner to help mitigate patient pressure and shear while in the lateral position. It is applicable for orthopedic, thoracic or urology procedures that dictate the patient be stable in the lateral position.

Vac-Pac® is a registered trademark of Natus Medical Incorporated

40612	Vac-Pac® Cover Medium Size 30	37 L x 28 W x ½" H	93.98 x 71.12 x 1.27 cm
40612-31	Vac-Pac® Cover Large Size 31	35 ½ L x 38 ¾ W x ½" H	90.17 x 98.43 x 1.27 cm
40613	Vac-Pac® Cover T-Shape Size 35	15 L x 34 W x ½" H (bottom) 9 L x 23 W x ½" H (top)	38.1 x 86.36 x 1.27 cm (bottom) 22.86 x 58.42 x 1.27 cm (top)