

Orthopedic Hip Positioner

- Frame-like positioning
- Developed at Rothman Institute



40218 Adult Wide Donut Head Pad 8 OD x 4 1/2 ID x 1 3/4" 20.5 OD x 11.4 ID x 4.5 cm



Peg Board Cover Set

- Covers hard peg board
- Small and Large post covers

40719 Peg Board and Post Pad Set Includes 40100, 4-40719A, 4-40719B

40719A Short Post Pad 5 x 5 x 1/2" 12.7 x 14 x 1.3 cm

40719B Tall Post Pad 5 x 8 1/2 x 1/2" 12.7 x 21.6 x 1.3 cm

Action® Head Pads

- 24 Head Pads to choose from
- Neonatal to Bariatric sizes
- Sample of styles below



40218 Adult Wide Donut Head Pad 8 OD x 4 1/2 ID x 1 3/4" 20.5 OD x 11.4 ID x 4.5 cm

40213 Contoured Head Pad 7 x 7 x 3 1/4" (1 1/2" contoured) 18 x 18 x 8 cm

40216L Large Prone Headrest 10 5/8 x 9 5 1/4" 27 x 22.0 x 13.4 cm

40204 Adult Horseshoe Head Pad 8 OD x 3 ID x 1 3/4" 20.5 OD x 7.5 ID x 4.5 cm

40207 Head Pad with Center Dish 9 x 10 x 1" 23 x 25 x 2.5 cm

Action® Arm and Lateral Positioners

- Arm Guards to protect skin integrity in Orthopedic cases
- For all positions

40130 Toboggan Pad 16 1/4 x 8 1/4 x 3/8" 41.3 x 21 x 1 cm

40303 Ulnar Forearm Protector 18 x 13 x 1/4" 46 x 33 x .6 cm

40302 Shoulder Brace 2 piece Set 5 x 6 x 1" 13 x 15 x 2.5 cm

40310 Medium Contoured Arm board Pad 20 x 5 3/4 x 2 1/4" 51 x 14.6 x 5.7 cm

40610 Lateral Positioner 28 x 19.5 x 6/1" 71.1 x 49.5 x 15.2/2.5 cm


Action® Specialty Frames

- McGuire, Relton, CHOP, Kambin, Wilson, Andrews and Montreal frame pads
- Maximize pressure and shear protection in vulnerable bony prominence areas
- Custom Manufacturing available
- Positioning Strap is fully lined in Akton® polymer- 2 sizes available




OPERATION DryMax®

- Single-use disposable product
- Safer work environment reduces slips, trips and falls
- Improves clean up efficiency for faster turnover times
- Improves infection control by trapping and retaining contaminated fluids
- Reduce red bag costs
- Lower laundry bills from cleaning towels and blankets used to soak up liquids
- Perfect fluid management for operating room, procedure or sink areas
- Great for Orthopedics, Urology, Labor & Delivery, Central Sterilization and more
- Adheres to OSHA, JCAHO and CMS standard regulations



DryMax 2.4 - superabsorbent

- Traps and retains fluids
- Low profile, ultra-thin mat
- Absorbs from top and bottom
- Safer work environment
- Fast, efficient turnover
- **9 liter capacity**




DryMax TRIPLE - superabsorbent

- Traps and retains fluids
- Flexible 3-part mat
- Full OR table length
- Low profile, ultra-thin mat
- Absorbs from top and bottom
- **15 liter capacity**











DryMax XL - absorbent

- Anti-skid design improves safety
- Large surface area keeps floors dry
- Can be placed under a stool
- Ideal to use under sink area
- Good for endoscopy based surgery and urology cases



DryMax COMBIMAT - flexible absorbent

- High absorbency combination of the 2.4 and XL
- Low profile, ultra-thin mat
- Anti-skid barrier
- Green mat absorbs from top and bottom
- Can easily be folded to protect shoes and floors
- Good for endoscopy based surgery, gynecology and urology cases

	Model	Size	Weight	Unique Quality	Packaging	Absorbency:  = 1 liter
DryMax 2.4	DM20241	28 x 14 1/2"	120 g	Absorbs from top and bottom	25 pcs/bag, 100/case	9 lt of water  3.5 lt of NaCl 
DryMax TRIPLE	DM20243	85 x 14 1/2"	200 g	Folds around objects	40/case	15 lt of water  6 lt of NaCl 
DryMax XL	DM22201	39 x 29 1/2"	160 g	Non-skid	10 pcs/bag, 50/case	1.5 lt of water  1.5 lt of NaCl 
DryMax COMBIMAT	DM22202	39 x 43"	350 g	2 mats in 1 XL + 2.4	5 pcs/bag, 30/case	15 lt of water  7 lt of NaCl 